

# END-TO-END POWER CHAIR SOLUTIONS



	JAZZY AIR® 2 <sup>1</sup>	JAZZY® CARBON <sup>1</sup>	JAZZY® PASSPORT <sup>1</sup>	GO CHAIR® 1 <sup>1</sup>	JAZZY® EVO 613 SERIES <sup>1</sup>	JAZZY® EVO 614 <sup>1</sup>	JAZZY® EVO 614HD <sup>1</sup>	GO CHAIR® MED <sup>1</sup>	JAZZY® ELITE ES <sup>1</sup>	JAZZY® SELECT <sup>1</sup>	JAZZY® ELITE 14 <sup>1</sup>	JAZZY® 600 ES <sup>1</sup>	JAZZY® ELITE HD <sup>1</sup>	JAZZY® 1450 <sup>1</sup>	
Weight Capacity	300 lbs.	300 lbs.	250 lbs.	300 lbs.	300 lbs.	300 lbs.	450 lbs.	300 lbs.	300 lbs.	300 lbs.	300 lbs.	300 lbs.	450 lbs. <sup>7</sup>	600 lbs.	
Turning Radius <sup>11</sup>	24.25"	24"	24.25"	25.5"	22"	22"	22"	25.5"	24.75"	22.5"	24"	20.5"	24"	26.5"	
Width <sup>11</sup>	23.25"	24"	23.5"	22.5"	22"	25.4"	25.4"	23"	22.75"	23.75"	24.25"	24.5"	24.25"	29.25"	
Length <sup>11,12</sup>	33.21"/42" (with footplate)	37" (with footplate)	40" (with footplate)	28.5"/33" (with footplate)	35" (with footplate)	40.5" (with footplate)	40.5" (with footplate)	28.5"/33" (with footplate)	34.75"/41" (with footplate)	32.99"/44" (with footplate)	35.5"/42.5" (with footplate)	35.5"/40.25" (with footplate)	35.5"	38.25"	
Maximum Speed <sup>3</sup>	Up to 4 mph (elevated and non-elevated)	Up to 3.7 mph	Up to 3.6 mph	Up to 3.8 mph	Up to 5.1 mph	Up to 4.4 mph	Up to 4.88 mph	Up to 3.7 mph	Up to 4 mph	Up to 4.3 mph	Up to 4 mph	Up to 4 mph	Up to 4 mph	Up to 4.9 mph	
Ground Clearance <sup>11</sup>	1" at motor bracket	2"	3" at motor	1.75" at motor	3" at motor bracket	3" at motor bracket	3" at motor bracket	1.75" at motor	1.85" at motors	.75" at motor bracket	3" at battery tray	2.75"	3"	1.5" at anti-tips	
Front Wheels	6" solid casters	7" solid casters	8" solid casters	6" solid casters	5" solid	6" solid	6" solid	6" solid casters	3" solid anti-tips	5" solid anti-tips	3" solid anti-tips	6" solid OMNI-casters	3" solid anti-tips	4" solid anti-tips	
Drive Wheels	10" solid	8.5" solid	11.25" perforated	9" solid	13" solid	13.75" perforated	13.75" perforated	9" solid	9" solid	10" solid	14" foam filled	14" foam filled	14" solid	14" solid	
Rear Wheels	6" solid casters	N/A	1.75" anti-tips	3" solid anti-tips	5" solid	6" solid	6" solid	3" solid anti-tips	6" solid casters	6" solid casters	6" solid casters	6" solid OMNI-casters	6" solid casters	9" solid casters	
Drivetrain	Two-motor, Mid-Wheel 6 <sup>®</sup> drive	24v/250W brushless motors	Two-motor, rear-wheel drive	Two-motor, rear-wheel drive	Two-motor, in-line, mid-wheel drive	Two in-line HD motors	Two in-line HD motors	Two-motor, rear-wheel drive	Two-motor, in-line, front-wheel drive	Two-motor, mid-wheel drive	Two-motor, front-wheel drive	Two-motor, mid-wheel drive	Two-motor, front-wheel drive	Two-motor, front-wheel drive	
Braking System	Regenerative and electro-mechanical	Regenerative and electro-mechanical	Regenerative and electro-mechanical	Regenerative and electro-mechanical	Regenerative and electro-mechanical	Regenerative and electro-mechanical	Regenerative and electro-mechanical	Regenerative and electro-mechanical	Regenerative and electro-mechanical	Regenerative and electro-mechanical	Regenerative and electro-mechanical	Regenerative and electro-mechanical	Regenerative and electro-mechanical	Regenerative and electro-mechanical	
Suspension Type	Active-Trac <sup>®</sup> (ATX) with Mid-Wheel 6 <sup>®</sup> Technology	Front	Limited	Limited	Active-Trac <sup>®</sup> (ATX)	Full independent front & rear Active-Trac <sup>®</sup> (ATX)	Full independent front & rear Active-Trac <sup>®</sup> (ATX)	Limited	Limited	Active-Trac <sup>®</sup> (ATX)	Limited	Active-Trac <sup>®</sup> (ATX) with Mid-wheel 6 <sup>®</sup> technology	Limited	Swing Arm	
High-Back Seat (max. size)	20"W x 20-22"D	N/A	N/A	N/A	20"W x 18-20"D	20"W x 18-20"D	18-22"W x 18-20"D	N/A	20"W x 20"D	20"W x 20"D	22"W x 22"D	20"W x 20"D	22"W x 22"D	32"W x 28"D	
High-Back Seat-to-Floor Range	21.25-33"	18.75-20.25" <sup>12</sup>	19"	N/A	19.75-21.25"	22-24"	22-24"	N/A	21.5-23.5"	21-22"	20-22"	20.5-22.5"	20-22"	20.5-22.5"	
Synergy <sup>®</sup> Seating	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	20"W x 20"D	20"W x 20"D	20"W x 20"D	24"W x 24"D	32"W x 28"D	
Synergy <sup>®</sup> Seat-to-Floor Range	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	17.25-19.25"	15.75-17.75"	15.875-17.875"	15.75-17.75"	17.25-19.25"	
Power Elevating Seat	Standard	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	
TRU-Balance <sup>®</sup>	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	
Specialty Controls	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	Yes	
Standard Electronics	40A, PG VR2	35A Brushless Micon	50A, PG, VSI	40A, PG GC3	40A, PG GC3	60A Dynamic Linx	60A Dynamic Linx	40A, PG GC3	40A, PG GC3	40A, PG GC3	50A, Dynamic LiNX	50A, Dynamic LiNX	50A, Dynamic LiNX	90A, PG VR2	
Battery Size	(2) 12 volt, U-1 (std.), 40 Ah (optional)	(1) 24V 12Ah (288WH) lithium ion	(1) 25.2V 18Ah lithium ion, (1) 25.2V 10.4 Ah lithium ion travel battery	(2) 18 Ah, 12 volt, deep-cycle	EVO 613 Li: (1) 25.6V, 32.4 Ah lithium iron phosphate EVO 613: (2) 12 volt, U-1	(2) NF-22	(2) NF-22	(2) 18 Ah, 12 volt, deep-cycle	(2) 12 volt, U-1	(2) 12 volt, U-1	(2) 12 volt, NF-22	(2) 12 volt, NF-22	(2) 12 volt, NF-22	(2) 12 volt, Group 24	
Standard Battery Charger	5A, off-board	2A, off-board	2A, off-board	2A, off-board	3.5A, off-board	5A, off-board	5A, off-board	2A, off-board	3.5A, off-board	3.5A, off-board	5A, off-board	5A, off-board	5A, off-board	8A, off-board	
Per-Charge Range <sup>3,6</sup>	Up to 15 miles at 300 lbs. Up to 16.6 miles at 220 lbs.	Up to 9.3 miles	Up to 5 miles (220 lbs.) (10.4 Ah) Up to 8.6 miles (220 lbs.) (18 Ah) Up to 4.6 miles (250 lbs.) (10.4 Ah) Up to 7.9 miles (250 lbs.) (18 Ah)	Up to 8.7 miles (300 lbs.) Up to 13.2 miles (200 lbs.)	EVO 613 Li: Up to 21 miles at 200 lbs. Up to 17.6 miles at 300 lbs. EVO 613: Up to 17.1 miles at 200 lbs. Up to 14 miles at 300 lbs.	Up to 19 miles at 300 lbs. Up to 21.5 miles at 200 lbs.	Up to 13 miles at 450 lbs. Up to 21.5 miles at 200 lbs.	Up to 8.6 miles	Up to 15 miles	Up to 19.4 miles	Up to 18 miles (300 lbs.) Up to 21 miles (200 lbs.)	Up to 24.8 miles (300 lbs.) Up to 29.5 miles (200 lbs.)	Up to 14.4 miles	Up to 16.5 miles	
Battery Weight <sup>8</sup>	24.5 lbs. (U-1), 32.5 lbs. each (40Ah)	4 lbs.	7.6 lbs. (18 Ah), 5.2 lbs (10.4 Ah travel battery)	14 lbs. each	EVO 613 Li: 18.5 lbs. EVO 613: 23.5 lbs. each (2)	36 lbs. each (2)	36 lbs. each (2)	14 lbs. each	25 lbs. each	25.2 lbs. each	37 lbs. each	37 lbs. each	37 lbs. each	52 lbs. each	
Weight of Heaviest Piece <sup>5</sup>	202 lbs. (with seat) <sup>9</sup>	39.6 lbs. <sup>10</sup>	59.8 lbs. (unit doesn't disassemble) <sup>10</sup>	36 lbs.	94 lbs. (base)	125.7 lbs. (base)	125.7 lbs. (base)	36 lbs.	69 lbs.	94 lbs.	99 lbs.	124 lbs.	99 lbs.	195 lbs.	
Standard Seat Weight	41 lbs. sport, high-back	N/A	N/A	30 lbs. (18") 32 lbs. (20") compact, foldable	46 lbs. (16") 48 lbs. (18") 50 lbs. (20") memory high-back	57.8 lbs.	57.8 lbs.	18" wide: 29 lbs. <sup>4</sup> 20" wide: 32 lbs. <sup>4</sup> compact, foldable	50.2 lbs. comfort, high-back	50.2 lbs. comfort, high-back	18" wide: 51.8 lbs. 20" wide: 54 lbs. 22" wide: 57.6 lbs. high-back	18" wide: 51.8 lbs. 20" wide: 54 lbs. high-back	18" wide: 51.8 lbs. 20" wide: 54 lbs. high-back	52 lbs. high-back	
Standard Seat Size	16" x 16-18", 18" x 18-20", 20" x 20-22"	17" x 17"	17.5" x 16.5"	18" x 17", 20" x 18"	16" x 16-18", 18" x 18-20", 20" x 18-20"	16" x 16-18", 18" x 18-20", 20" x 18-20"	18" x 18-20", 20" x 18-20", 22" x 20-22"	18" x 17", 20" x 18"	16" x 16-18", 18" x 18-20", 20" x 18-20"	16" x 16-18", 18" x 18-20", 20" x 18-20"	18" x 18-20", 20" x 18-20", 22" x 20-22"	18" x 18-20", 20" x 18-20"	18" x 18-20", 20" x 18-20", 22" x 20-22", 24" x 24"	24" x 24"	
HCPCS Codes	K0899 POWER MOBILITY DEVICE NOT CODED BY DME PDAC OR DOES NOT MEET CRITERIA	K0899 POWER MOBILITY DEVICE NOT CODED BY DME PDAC OR DOES NOT MEET CRITERIA	K0899 POWER MOBILITY DEVICE NOT CODED BY DME PDAC OR DOES NOT MEET CRITERIA	K0899 POWER MOBILITY DEVICE NOT CODED BY DME PDAC OR DOES NOT MEET CRITERIA	K0822/K0823 Group 2: Standard	K0822/K0823 Group 2: Standard	K0824/K0825 Group 2: Heavy Duty	K0821 Group 2: Standard, Portable	K0822/K0823 Group 2: Standard	K0824/K0825 Group 2: Heavy Duty	K0826/K0827 Group 2: Very Heavy Duty K0839 Single Power Very Heavy Duty				
Warranty	5-year limited on frame; 13-month limited on electronics; 13-month on drive motors; 13-month on batteries	Lifetime limited on frame; 1-year limited on electrical control/lever system, motor/drive system, bearing sleeve and battery	5-year limited on frame; 13-month limited on electronics; 13-month on drive motors; 13-month on batteries	5-year limited on frame; 13-month limited on electronics; 13-month on drive motors; 13-month on batteries	5-year limited on frame; 13-month limited on electronics; 13-month limited on SLA batteries; 3-year limited on lithium batteries	5-year limited on frame; 13-month limited on electronics; 13-month limited on drive motors; 13-month limited on batteries	5-year limited on frame; 13-month limited on electronics; 13-month limited on drive motors; 13-month limited on batteries	5-year limited on frame; 13-month limited on electronics; 13-month on drive motors; 13-month on batteries	5-year limited on frame; 13-month limited on electronics; 13-month on drive motors; 13-month on batteries	5-year limited on frame; 13-month limited on electronics; 13-month on drive motors; 13-month on batteries	5-year limited on frame; 13-month limited on electronics; 13-month on drive motors; 13-month on batteries	5-year limited on frame; 13-month limited on electronics; 13-month on drive motors; 13-month on batteries	Lifetime limited on frame; 13-month limited on electronics; 13-month on drive motors; 13-month on batteries	5-year limited on frame; 13-month limited on electronics; 13-month on drive motors; 13-month on batteries	K0826/K0827 Models: Lifetime limited on frame; 13-month limited on electronics; 13-month limited on drive motors; 13-month on batteries K0839 Model: Lifetime limited on frame; 2-year limited on electronics; 18-month limited on drive motors; 13-month on batteries
<b>NOT CODED BY DME PDAC</b>															

1) This product is a FDA Class II Medical Device designed to aid individuals with mobility impairments.  
 2) Without front riggings.  
 3) Varies with user weight, terrain, battery type, battery condition and tire pressure.  
 4) Seat weight includes joystick.  
 5) Does not include seat or battery weights.  
 6) Tested in accordance with ANSI/RESNA, WC Vol. 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Testing conducted at maximum weight capacities. This specification can be subject to a variance of +10%, -5%.

7) Maximum weight capacity is reduced to 400 lbs. (181 kg) on 18" deep seat configurations.  
 8) Battery weight may vary based on manufacturer.  
 9) Seat weight included (because tools are required to remove seat). Battery weight not included.  
 10) Seat weight included. Unit does not disassemble.  
 11) Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.  
 12) With and without standard cushion.

BATTERY HCPCS CODES	
U-1	E2365
NF-22	E2361
Group 24	E2363
Lithium	E2397
18 AH	K0733

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. The information contained herein is correct at time of publication; we reserve the right to alter specifications without prior notice.

