



MODEL	LC-358PW ^{1,4}	LC-358S ^{1,4}	LC-358M ^{1,4}	LC-358L ^{1,4}	LC-358XL ^{1,3,4}	LC-358XXL ^{1,3,5}	LC-105 ¹	LC-250 ¹
Size	Petite Wide	Small	Medium	Large	Heavy-Duty	Very Heavy-Duty	Small	Medium
Weight Capacity	375 lbs.	375 lbs.	375 lbs.	375 lbs.	500 lbs.	600 lbs.	325 lbs.	375 lbs.
Suggested Height Range ⁶	5'3" & Below	5'3" & Below	5'4" - 5'9"	5'9" - 6'	5'8" - 6'1"	5'8" - 6'1"	5'3" & Below	5'4" - 5'9"
Position	3-Position	3-Position	3-Position	3-Position	3-Position	2-Position	3-Position	3-Position
Overall Width	40"	35"	35"	37"	39.5"	43"	33"	34"
Chair Weight	130 lbs.	118 lbs.	117 lbs.	120 lbs.	159 lbs.	168 lbs.	96 lbs.	118 lbs.
Seating Type	Total Comfort ²	Total Comfort ²	Total Comfort ²	Total Comfort ²	Sinuous Spring	Sinuous Spring	Sinuous Spring	Total Comfort ²
Back Style	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Sewn Pillow Back	Blown Fiber Split Back
Seat Width	24"	18"	20"	22"	26"	30"	20"	22"
Seat Depth	19"	19"	20"	22"	22"	22"	19"	20"
Seat to Floor Height	19"	19"	20"	21"	20"	20"	18"	20"
Top of Back to Seat	23.5"	24"	24.5"	25.5"	29"	29"	24"	26"
Distance From Wall	18"	18"	17"	19"	19"	14"	17"	20"
Fully Padded Chaise	Yes	Yes	Yes	Yes	Yes (Removable)	Yes (Removable)	Yes	Yes
Footrest Extension	Standard	N/A	Standard	Standard	N/A	N/A	N/A	N/A
Additional Options	Heat & Massage ⁴ , Lithium Battery Backup	Heat & Massage ⁴ , Lithium Battery Backup	Heat & Massage ⁴ , Lithium Battery Backup	Heat & Massage ⁴ , Lithium Battery Backup	Heat & Massage ⁴ , Lithium Battery Backup	Lithium Battery Backup	No	No
Warranty	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	1-Yr. Limited	2-Yr. Limited

FABRIC SELECTION GUIDE	Standard	Optional
	<p>Cloud 9^{QS}</p> <p>BLACK CHERRY, PACIFIC, STONE, WALNUT</p> <p>Cloud 9^{QS}</p> <p>STONE, WALNUT</p> <p>Micro-Suede^{QS}</p> <p>COCOA, SANDAL</p> <p>Cloud 9^{QS}</p> <p>PACIFIC, STONE, WALNUT</p>	<p>Crypton Aria[®]</p> <p>ESPRESSO, SAND, COOL GRAY, RED</p> <p>Lexis Sta-Kleen[®]</p> <p>BUFF, FUDGE, GARNET, PECAN, CHARCOAL, CHESTNUT, BLACK, BURGUNDY, MUSHROOM</p> <p>Lexis Sta-Kleen[®]</p> <p>CHESTNUT</p> <p>Lexis Urethane[®]</p> <p>NEW CHESTNUT</p> <p>Lexis Sta-Kleen[®]</p> <p>CHESTNUT</p>

• Recommended HCPCS: E0627—Seat lift mechanism, electric, any type. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

• Made of foreign and domestic components

FOOTNOTES:

- Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- Unique spring design and a one-piece soft foam seat for outstanding relaxation.
- Features Dual Motor Sync System.
- USB hand control not available with Heat & Massage options.
- Only available in Cloud 9 Stone, Cloud 9 Walnut, and Lexis Sta-Kleen® Chestnut.
- Height ranges should be used only as a guide due to varying factors of each individual user.

Like all furniture, Pride® Power Lift Recliners are individually built and there may be differences in actual measurements due to variations in fabric, upholstery, and filling, in addition to the method of measurement.

ADDITIONAL NOTES:

- All chairs feature premium Furniture-Grade Laminated/Hardwood frames.
- All chairs have a Lifetime Warranty on lift mechanism, shell and recline mechanism (except LC-358XL, LC-358XXL, & LC-105).
- All chairs have integrated emergency battery backup (batteries not included).
- All chairs can have fiber fill added or removed for custom comfort adjustment.
- All chairs are shipped in an impact-resistant, reinforced, corrugated box.

POSITION GUIDE	
PRODUCT VIEW	
2-Position	The 2-position TV recline is achievable on all Pride® Power Lift Recliners. This is a comfortable position for performing many relaxing activities including reading, watching TV, working on your tablet or for conversation with family and friends.
3-Position	The 3-position full recline is available on most Pride® Power Lift Recliner models and puts you in a position ideal for napping or relaxing.

The lifted, seated and reading positions are achievable on all Pride® Power Lift Recliners.

Rev B/February 2024

I NFI NFB4553

NOTE: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

PARTNER TRADEMARKS



ICON KEY

