

**FOOTNOTES**

1. Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
2. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.
3. Battery weight not included.
4. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%.
5. Battery pack is 19.5 lbs. with 12AH batteries.

6. Deep cycle (AGM or gel cell type recommended).
7. 18 AH battery pack weighs 29.5 lbs.
8. Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
9. Battery weight may vary based on manufacturer.
10. Width of base (not including seat).
11. Includes standard seat weight.

12. Must be billed to Medicare as either a beneficiary requested upgrade with a properly executed ABN, or as an upgrade at no additional charge to the beneficiary.
13. Light package.
14. Interchangeable shrouds (all shroud panels included).
15. Results derived from theoretical calculation based on battery specifications and drive system performance. Test conducted at maximum weight capacity.



END-TO-END



	iRIDE® 2 <small>Non-Medical Device*</small>	i-Go™ <small>Non-Medical Device*</small>	Go Go® Endurance Li <small>FDA Class II Medical Device¹</small>	Go Go® Ultra X <small>FDA Class II Medical Device¹</small>	Go Go® ES 2 <small>FDA Class II Medical Device¹</small>	Go Go® Elite Traveller® <small>FDA Class II Medical Device¹</small>	Go Go® Sport <small>FDA Class II Medical Device¹</small>	Go Go® LX w/ CTS Suspension <small>FDA Class II Medical Device¹</small>	Revo® 2.0 <small>FDA Class II Medical Device¹</small>	Zero Turn 8 <small>FDA Class II Medical Device¹</small>	Zero Turn 10 <small>FDA Class II Medical Device¹</small>	Pursuit 2 <small>FDA Class II Medical Device¹</small>	PX4 <small>FDA Class II Medical Device¹</small>	Victory® 9 <small>FDA Class II Medical Device¹</small>	Victory® 10 <small>FDA Class II Medical Device¹</small>	Victory® 10.2 <small>FDA Class II Medical Device¹</small>	Victory® LX Sport¹³ <small>FDA Class II Medical Device¹</small>	Maxima¹³ <small>FDA Class II Medical Device¹</small>	
<b>MODEL NUMBER</b>	S25	SC20	554XLIT	S39 (3-Wheel) S49 (4-Wheel)	581	SC40E (3-Wheel) SC44E (4-Wheel)	S73 (3-Wheel) S74 (4-Wheel)	S50LX	S66 (3-Wheel) S67 (4-Wheel)	JAZZYZTS	S710ZT	SC7132	SC134	SC609 (3-Wheel) SC709 (4-Wheel)	SC610 (3-Wheel) SC710 (4-Wheel)	S6102	S710LXW	SC901 (3-Wheel) SC941 (4-Wheel)	
<b>WEIGHT CAPACITY</b>	250 lbs.	275 lbs.	325 lbs.	300 lbs.	250 lbs.	300 lbs.	325 lbs.	300 lbs.	400 lbs.	325 lbs.	400 lbs.	400 lbs.	500 lbs.	300 lbs.	400 lbs.	400 lbs.	400 lbs.	400 lbs.	500 lbs.
<b>SPEED (MAXIMUM)²</b>	Up to 5.8 mph	Up to 4 mph	Up to 4.1 mph	3-Wheel: Up to 4.3 mph 4-Wheel: Up to 4.1 mph	Up to 4 mph	Up to 4 mph	Up to 4.7 mph	Up to 6 mph	Up to 5.2 mph	Up to 5.9 mph	Up to 7.3 mph	Up to 9.3 mph	Up to 9.2 mph	Up to 5.3 mph	Up to 5.3 mph	Up to 5.2 mph	Up to 8 mph	Up to 5.8 mph	
<b>GROUND CLEARANCE⁴</b>	1.3" at motor	4.2" at mid-frame	1.5" at motor 3.25" at side deck	1.25" at motor 3.2" at mid-frame	0.75" at anti-tip 1.5" at motor	0.75" at anti-tip 2.25" at mid-frame	1.5" at motor 3" at mid-frame	0.75" at anti-tip 1.5" at motor 2.5" at frame	1.25" at transaxle 3.25" at mid-frame (3-Wheel) 3.5" at mid-frame (4-Wheel)	1.3" at motor 1.5" at anti-tip 3.5" at side deck	2" at motor 3" at mid frame	3.78" at transaxle 4.6" at mid-frame	3.5" at motor 4.4" at mid-frame	3-Wheel: 1.4" at anti-tip 1.8" at motor 2" at mid-frame 4-Wheel: 1.25" at anti-tip 1.75" at motor 2.5" at mid-frame	3-Wheel: 1.26" at anti-tip 2.2" at motor 3" at mid-frame 4-Wheel: 2" at motor 1.5" at anti-tip 3" at mid-frame	2.4" at transaxle 2.1" at anti-tip 3" at side deck	2.2" at motor 3.15" at rear mid-frame	3-Wheel: 2.5" at motor 3" at mid-frame 4-Wheel: 2.25" at motor	
<b>3-WHEEL TURNING RADIUS⁴</b>	33.5"	45.9"	N/A	36.5"	31"	33"	36"	33.5"	40"	N/A	N/A	N/A	N/A	40.25"	42"	42"	N/A	49"	
<b>4-WHEEL TURNING RADIUS⁴</b>	N/A	N/A	48.25"	52.25"	N/A	44"	50.5"	N/A	52"	37.25"	43"	69.3"	82.3"	51.5"	54"	N/A	61.3"	53"	
<b>3-WHEEL LENGTH/WIDTH⁴¹⁰</b>	38"/19.25"	44"/20.3"	N/A	40.25"/19.25"	36.75"/19.5"	37"/19.5"	42"/21.25"	39.75"/20.5"	46"/20.75"	N/A	N/A	N/A	N/A	43.5"/22.25"	45.25"/22.25"	45"/22.25"	N/A	47.5"/27"	
<b>4-WHEEL LENGTH/WIDTH⁴¹⁰</b>	N/A	N/A	42"/20.5"	40.25"/19.25"	N/A	39.5"/19.5"	43"/21.25"	41.75"/20.5"	47"/20.75"	45.75"/20.75"	48"/24.25"	57.5"/29.7"	56"/27" - 27.8"	45.5"/22.5"	47"/22.25"	N/A	47"/25.5"	47.5"/27"	
<b>FRONT TIRES⁴</b>	6.5" solid	8" hub motor	8" solid	8" solid	7.5" solid	3-Wheel: 8" solid 4-Wheel: 9" solid	3-Wheel: 8" solid 4-Wheel: 9" solid	8" solid	3-Wheel: 8" solid 4-Wheel: 9" solid	7" solid	9" solid	13" pneumatic	13" pneumatic	3-Wheel: 8" solid 4-Wheel: 9" solid	10" solid	10" solid	10" solid	10" solid	3-Wheel: 10.5" solid 4-Wheel: 10" solid
<b>REAR TIRES⁴</b>	6.75" solid	9" perforated	8" solid	8" solid	7.5" solid	8" solid	9" solid	8" solid	9" solid	9" solid	10.75" solid	13" pneumatic	13" pneumatic	9" solid	10" solid	10" solid	10" solid	10" solid	10.75" solid
<b>SUSPENSION</b>	N/A	N/A	Front & Rear	N/A	N/A	N/A	N/A	Front & Rear	Front & Rear	Front & Rear	Front & Rear	Front & Rear	Front & Rear	N/A	N/A	N/A	Front & Rear	N/A	
<b>RANGE PER CHARGE²,⁸</b>	Up to 12 miles (200 lbs.) Up to 10 miles (250 lbs.)	Up to 9.3 miles (220 lbs.) Up to 7 miles (275 lbs.)	Up to 6.6 miles (200 lbs, 8 AH) Up to 13.2 miles (200 lbs, 16 AH)	3-Wheel: Up to 8 miles (200 lbs.) Up to 6.2 miles (300 lbs.) 4-Wheel: Up to 8.9 miles (200 lbs.) Up to 6.6 miles (300 lbs.)	Up to 6.24 miles (250 lbs.)	3-Wheel: Up to 7.3 miles (12Ah, 300 lbs.) Up to 10.5 miles (18Ah, 300 lbs.) 4-Wheel: Up to 6 miles (12Ah, 300 lbs.) Up to 9.75 miles (18Ah, 300 lbs.)	3-Wheel: Up to 12.75 miles (325 lbs.) 4-Wheel: Up to 12 miles (325 lbs.)	Up to 6 miles (300 lbs.)	3-Wheel: Up to 17.4 miles (200 lbs.) Up to 12.8 miles (400 lbs.) 4-Wheel: Up to 17.8 miles (200 lbs.) Up to 14.8 miles (400 lbs.)	Up to 13.1 miles (200 lbs.) Up to 11 miles (325 lbs.)	Up to 21.2 miles (200 lbs.) Up to 15.9 miles (400 lbs.)	Single Battery: Up to 17.4 miles (220 lbs.) Up to 13.8 miles (400 lbs.) Dual Battery: Up to 34.8 miles (220 lbs.) Up to 30.9 miles (350 lbs.)	Up to 23.4 miles (220 lbs.) Up to 14 miles at (500 lbs.)	3-Wheel: Up to 18 miles (300 lbs.) 4-Wheel: Up to 16.3 miles (300 lbs.)	Up to 15.5 miles (400 lbs.)	Up to 13 miles (40Ah, 400 lbs.) Up to 11 miles (J-1, 400 lbs.)	Up to 19.7 miles (200 lbs.) Up to 13 miles (400 lbs.)	3-Wheel: Up to 19.8 miles (500 lbs.) 4-Wheel: Up to 19.5 miles (500 lbs.)	
<b>TOTAL WEIGHT WITHOUT BATTERIES³,⁴,¹¹</b>	42 lbs.	60.6 lbs.	93 lbs.	3-Wheel: 76.3 lbs. 4-Wheel: 81.5 lbs.	58 lbs.	3-Wheel: 75.5 lbs. 4-Wheel: 79.5 lbs.	3-Wheel: 92.75 lbs. 4-Wheel: 100 lbs.	3-Wheel: 85 lbs. 4-Wheel: 90 lbs.	3-Wheel: 120 lbs. 4-Wheel: 132 lbs.	120 lbs.	181.8 lbs.	220.6 lbs.	253.8 lbs.	3-Wheel: 117 lbs. 4-Wheel: 131 lbs.	3-Wheel: 133 lbs. 4-Wheel: 137 lbs.	130.25 lbs.	160 lbs.	3-Wheel: 183.5 lbs. 4-Wheel: 196.5 lbs.	
<b>WEIGHT OF HEAVIEST PIECE⁴</b>	32 lbs. (Scooter Base)	32 lbs. (Rear Base with Seat Back)	40.75 lbs. (Front Section)	3-Wheel: 29.5 lbs. (Front Section) 4-Wheel: 34.5 lbs. (Front Section)	41.4 lbs. (Scooter Base)	3-Wheel: 29.5 lbs. (Front Section) 4-Wheel: 34 lbs. (Front Section)	3-Wheel: 35.75 lbs. (Front Section) 4-Wheel: 44 lbs. (Front Section)	33 lbs. (Front Section)	3-Wheel: 48 lbs. (Rear Section) 4-Wheel: 57 lbs. (Front Section)	48 lbs. (Rear Section)	75 lbs. (Front Section)	163 lbs. (Base)	201.2 lbs. (Base)	3-Wheel: 48.5 lbs. (Front Section) 4-Wheel: 59 lbs. (Front Section)	3-Wheel: 54 lbs. (Front Section) 4-Wheel: 63 lbs. (Front Section)	56.5 lbs. (Front Section)	77 lbs. (Front Section)	N/A (unit does not disassemble)	
<b>STANDARD SEATING</b>	Compact, Foldable (14.75" x 13.5")	Compact, Foldable with Removable Back (16" x 16")	Compact, Foldable (17" x 17")	Compact, Mesh, Foldable (17" x 17")	Compact, Foldable (17" x 15")	Compact, Foldable (17" x 17")	Compact, Foldable (17" x 17")	Compact, Foldable (17" x 17")	Compact, Foldable (20" x 18")	Compact, Foldable (17" x 17")	Compact, Foldable (18" x 17")	Limited Recline High Back, with Sliders (18" X 18")	Limited Recline High Back (18"x18")	Compact, Foldable (18" x 17")	Compact, Foldable (18" x 17")	Foldable (18"x17")	Limited Recline High Back with Sliders & Headrest (18" x 18")	Depth-Adjustable Limited Recline High Back (22" x 18"-20")	
<b>BATTERY REQUIREMENTS</b>	(1) 36V, 5.2Ah Lithium-ion	(1) 25.2V, 9.6Ah Lithium-ion	8Ah Lithium-ion pack or 16Ah Lithium-ion Pack	(2) 12Ah (SLA/AGM)	(2) 12Ah (SLA/AGM)	(2)12Ah (SLA/AGM) or (2)18Ah (SLA/AGM)	(2) 18Ah (SLA/AGM)	(2) 12Ah (SLA/AGM)	(2) U-1 (SLA/AGM)	(2) 22Ah (SLA/AGM)	(2) 40Ah (SLA/AGM)	(1) or (2) 24V, 50Ah Lithium-ion	(2) 12V 75Ah (SLA/AGM)	(2) U-1 (SLA/AGM)	(2) U-1 (SLA/AGM) or (2) 40Ah (SLA/AGM)	(2) U-1 (SLA/AGM) or (2) 40Ah (SLA/AGM)	(2) 50Ah (SLA/AGM)	(2) GEL, 58Ah, NF-22 (Nano-silica)	
<b>BATTERY WEIGHT⁹</b>	4.5 lbs.	4.8 lbs.	8 lbs. (8Ah)/ 12 lbs. (16Ah)	9 lbs. ea.⁵	9 lbs. ea.⁵	9 lbs. ea.⁵ / 14 lbs. ea.⁷	14 lbs. ea.⁷	9 lbs. ea.⁵	25 lbs. ea.	32 lbs. (Battery box)	30.25 lbs. ea.	23.8 lbs. ea.	48.8 lbs. ea.	24 lbs. ea.	24 lbs. ea./ 32 lbs. ea.	24 lbs. ea./ 32 lbs. ea.	30.2 lbs. ea.	36.5 lbs. ea.	
<b>BATTERY CHARGER</b>	Off-board, 36V (Lithium-ion)	Off-board, 3.5A (Lithium-ion)	Off-board, 3.5A (Lithium-ion)	Off-board, 2A	Off-board, 2A	Off-board, 2A	Off-board, 2A	Off-board, 2A	Off-board, 3.5A	Off-board, 3.5A	Off-board, 3.5A	Off-board, 8A	Off-board, 8A	Off-board, 3.5A	Off-board, 3.5A (U-1)/ Off-board, 5A (40Ah)	Off-board, 3.5A (U-1)/ Off-board, 5A (40Ah)	Off-board, 5A	Off-board, 5A	
<b>SEAT-TO-DECK HEIGHT RANGE</b>	11.5" - 12.5"	14.6"	16.25"-17.75"	16.5"-18"	16.25"-17.25"	15.75"-17.25"	17"-18"	16.25"-17.75"	16.5" - 18"	17"-18.5"	17.75"-19.25"	17.9"-20.9"	16.54"-18.39"	17.25"-18.75"	17.5"-19"	17.5" - 18"	16.5"-18.5"	17.13" - 19.13"; 17.36"-23.61" (with power elevating seat)	
<b>SEAT-TO-GROUND HEIGHT RANGE⁴</b>	19" - 20"	22.25"	21.25"-22.75"	21.5" - 23"	20" - 21"	18.25"-20.25"	20.5"-22"	20.75"-22.25"	22" - 23.5"	22" - 23.5"	23.25"-24.75"	28.74"	24.17"-25.59"	22.5" - 24"	23.5" - 25"	23.5" - 25"	23" - 25"	22.1" - 24.1" 22.32"-28.57" (with power elevating seat)	
<b>STANDARD BASE COLORS</b>	Mango Arctic Ice Blackberry Raspberry	Sugar Plum (matte) Iceberg Blue (matte) Black (matte) White (matte)	Garnet Red (Matte) Euro Grey¹⁴	Black	Black	Red Blue¹⁴	Red Blue¹⁴	Red Blue¹⁴	Red Blue¹⁴	Grey Street, True Blue	Green Moonstone Onyx Black Pearl White Sapphire Blue Fire Opal, Pewter	Candy Apple Red Ocean Blue (matte) Pearl White Black (matte)	Black Grey (matte) Orange	Satin Aluminum Dark Violet Red Peacock Blue Sunflower	Candy Apple Red	Candy Apple Red Viper Blue	Candy Apple Red Ocean Blue (matte)	Candy Apple Red True Blue	Garnet Red (matte) Ocean Blue (matte)
<b>HCPCS CODE</b>	Non-Medical*	Non-Medical*	K0899**	K0899**	K0800 - Group 1 Standard (3-Wheel)	K0800 - Group 1 Standard	K0801 - Group 1 Heavy Duty	K0899**	K0899**	K0899**	K0899**	K0899**	K0899**	K0800 - Group 1 standard	K0801 - Group 1 Heavy Duty	K0801 - Group 1 Heavy Duty	K0899**	K0802- Group 1 Very Heavy Duty (3-Wheel) K0808- Group 2 Very Heavy Duty (4-Wheel)¹²	
<b>WARRANTY</b>	3-year limited on frame 2-year limited on electronics 1-year limited on Drivetrain, Motor, Battery, and Charger	3-year limited on Frame 2-year limited on Electronics 1-year limited on Drivetrain, Motor, Battery, and Charger	Lifetime limited on Frame; 2-year limited on Electronics 2-year limited on Drivetrain 2-year limited on Battery	Lifetime limited on Frame 2-year limited on Electronics 2-year limited on Drivetrain 6-month limited on Batteries	Lifetime limited on Frame 2-year limited on Electronics 2-year limited on Drivetrain 6-month limited on Batteries	Lifetime limited on Frame 2-year limited on Electronics 2-year limited on Drivetrain 6-month limited on Batteries	Lifetime limited on Frame 2-year limited on Electronics 2-year limited on Drivetrain 6-month limited on Batteries	Lifetime limited on Frame 2-year limited on Electronics 2-year limited on Drivetrain 6-month limited on Batteries	Lifetime Limited on Frame; 2-year limited on Electronics; 2-year limited on Drivetrain 6-month limited on Batteries	Lifetime limited on Frame 2-year limited on Electronics 2-year limited on Drivetrain 6-month limited on Batteries	Lifetime limited on Frame 2-year limited on Electronics 2-year limited on Drivetrain 6-month limited on Batteries	Lifetime Limited on Frame 2-year limited on Electronics 2-year limited on Drivetrain 6-month limited on Batteries	Lifetime limited on Frame 2-year limited on Electronics 2-year limited on Drivetrain 3-year limited on Battery	Lifetime limited on Frame 2-year limited on Electronics 2-year limited on Drivetrain 6-month limited on Battery	Lifetime Limited on Frame 2-year limited on Electronics 2-year limited on Drivetrain 6-month limited on Batteries	Lifetime limited on Frame 2-year limited on Electronics 2-year limited on Drivetrain 6-month limited on Batteries	Lifetime limited on Frame 2-year limited on Electronics 2-year limited on Drivetrain 6-month limited on Batteries	Lifetime limited on Frame 2-year limited on Electronics 2-year limited on Drivetrain 6-month limited on Batteries	Lifetime limited on Frame 2-year limited on Electronics 2-year limited on Drivetrain 6-month limited on Batteries

**TRAVEL MOBILITY**

**ZERO-TURN SERIES**

**MID-SIZE LUXURY**

**HEAVY-DUTY**

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%.

DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage, and bundling guidelines may apply.

\*Non-Medical Devices: This product is not a medical device and is not intended to assist, treat, diagnose, or alleviate any medical condition or disability.  
\*\*K0899: Power mobility device, not coded by DME PDAC.