

THE MOST COMFORTABLE RECLINER IN THE WORLD!



COLLECTION	ULTRA ^{1,2,3}	ULTRA ^{1,2,3}	ULTRA ^{1,2,3}	RADIANCE ^{1,2,3}	RADIANCE ^{1,2,3}	RADIANCE ^{1,2,3}	RADIANCE ^{1,2,3}	ATLAS PLUS 2 ^{1,2}	TRANQUIL 2 ¹	TRANQUIL 2 ¹	TRANQUIL 2 ¹	ELEGANCE 2 ¹	ELEGANCE 2 ¹	METRO 2 ¹	URBANA 2 ¹
MODEL	PLR-4955S	PLR-4955M	PLR-4955L	PLR-3955PW	PLR-3955S	PLR-3955M	PLR-3955LT	PLR-2985M	PLR-935S	PLR-935M	PLR-935LT	PLR-975M	PLR-975L	PLR-925M	PLR-965M
SIZE	Small	Medium	Large	Petite Wide	Small	Medium	Large/Tall	Medium	Small	Medium	Large/Tall	Medium	Large	Medium	Medium
WEIGHT CAPACITY	400 lbs.	400 lbs.	400 lbs.	400 lbs.	400 lbs.	400 lbs.	400 lbs.	400 lbs.	400 lbs.	400 lbs.	400 lbs.	400 lbs.	400 lbs.	400 lbs.	400 lbs.
POSITION**	True-Infinite 5'4" and below	True-Infinite 5'5" - 6'	True-Infinite 6' and above	True-Infinite 5'4" and below	True-Infinite 5'4" and below	True-Infinite 5'5" - 6'1"	True-Infinite 6'2" and above	True-Infinite 5'4"-5'9"	True-Infinite 5'3" and below	True-Infinite 5'4"-5'9"	True-Infinite 5'10" and above	True-Infinite 5'4"-5'9"	True-Infinite 5'10" and above	True-Infinite 5'3" and below	True-Infinite 5'3" and below
OVERALL WIDTH	37"	39"	39"	41"	37"	39"	40"	39"	33.75"	37.5"	38"	35.5"	37.5"	33"	33"
CHAIR WEIGHT	158 lbs.	158 lbs.	176 lbs.	161 lbs.	156 lbs.	158 lbs.	178 lbs.	164 lbs.	137 lbs.	145 lbs.	155 lbs.	152 lbs.	157 lbs.	142 lbs.	145 lbs.
CHAIR HEIGHT	41"	43"	47"	40"	41"	43"	49"	42"	40"	42"	47"	42"	43"	41"	42"
SEATING TYPE	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort
SEAT WIDTH	20"	22.5"	23"	24"	20"	22.5"	23.5"	22"	20"	22"	23.5"	21"	22"	21"	21"
SEAT DEPTH	18"	21"	22"	18"	18"	21"	22"	21"	18.5"	20.5"	22.5"	21"	22"	21"	20"
SEAT TO FLOOR HEIGHT	18"	19.5"	21.5"	19.5"	18"	19.5"	21.5"	19.5"	18"	19.5"	21.5"	20"	20.5"	19"	20"
TOP OF BACK TO SEAT	25"	28"	30"	23"	25"	28"	31"	26"	25"	27"	30"	27"	29"	26.5"	25"
DISTANCE FROM WALL	15"	15"	21.5"	14"	15"	15"	21.5"	16"	13"	16"	16"	14"	16"	14"	14"
FULLY PADDED CHAISE	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
FOOTREST EXTENSION	Standard	Standard	Standard	Standard	Standard	Standard	Standard	Standard	Standard	Standard	Standard	Standard	Standard	Standard	Standard
HEAT	Standard	Standard	Standard	Standard	Standard	Standard	Standard	NA	NA	NA	NA	NA	NA	NA	NA
WARRANTY	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated

<p>STANDARD FABRIC SELECTION GUIDE</p>	<p>CAPRICCIO</p> SLATE	<p>CANYON</p> OCEAN SILT STEEL WALNUT	<p>BADLANDS</p> MUSHROOM WALNUT STEEL	<p>ASTRO</p> GREY BROWN MUSHROOM	<p>BADLANDS</p> MUSHROOM WALNUT STEEL	<p>SAVILLE</p> GREY BROWN	<p>STONEWASH</p> GRANITE GUNMETAL

STANDARD FEATURES:

- POWER HEADREST
- POWER LUMBAR
- USB REMOTE
- LITHIUM BATTERY BACKUP
- FOOTREST EXTENSION

POSITION GUIDE*:

- RELAXED
- RECLINED
- UPWARD LIFT
- TRUE INFINITE

= Indicates shortest available lead time.

1. Pride[®] FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
 2. The wireless charger is only compatible with smart phones that have wireless charging capabilities. Additionally, no other smart devices, (i.e. smart watches) should be used with the wireless charger.
 3. Unit has three different heating pads in the shoulders, lumbar and seat.

*Not all units achieve all positions. Please refer to specific product.
 **Height ranges should be used as a guide only due to varying factors of each individual user.
 Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.