

Model	Go Go Elite Traveller® Plus ¹
Part Number	3-wheel: SC53 4-wheel: SC54
Weight Capacity	300 lbs.
Maximum Speed ²	Up to 4.5 mph
Ground Clearance ⁴	3-wheel: 1.7" at transaxle, 0.83" at anti-tip, 2.6" at side deck
Turning Radius ⁴	3-wheel: 34.5" 4-wheel: 45"
Overall Length ⁴	40"
Overall Width ⁴	21.25"
Seat-to-Ground Height	20.5" - 22"
Seat-to-Deck Height ³	16.25"-17.75"
Front Tires	2" x 8" solid
Rear Tires	3" x 9" solid
Suspension	N/A
Range Per Charge ^{2,3}	3-Wheel: Up to 10.25 miles 4-Wheel: Up to 14.5 miles
Total Weight With Batteries	3-wheel: 116 lbs. 4-wheel: 116.5 lbs.
Total Weight Without Batteries	3-wheel: 86.5 lbs. 4-wheel: 87 lbs.
Heaviest Piece When Disassembled ⁴	3-wheel: 32.5 lbs. (rear section) 4-wheel: 35 lbs. (front section)
Standard Seating	Type: Foldable, compact, vinyl-covered Weight: 23 lbs. Material: Black vinyl Dimensions: Width: 17" Depth: 17"
Motor Output	24-volt DC motor, sealed mini-transaxle
Brakes	Regenerative and electromechanical
Battery Requirements ⁵	Size: (2) 18 AH Weight: 14 lbs. each, 29.5 lbs. per battery pack
Battery Charger	Off-board, 2 amp

HCPCS Code: K0806 (Group 2, Standard)

- 1 Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- 2 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- 3 Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 4 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- 5 Battery weight may vary based on manufacturer.

(US) 800-800-8586 • (Canada) 888-570-1113

pridemobility.com |    

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

0004952-102622

GO GO
ELITE TRAVELLER® PLUS

**FEATURES...
PERFORMANCE...
...AND
TRAVEL READY!**



FDA Class II Medical Device¹

LIVE YOUR BEST®
pridemobility.com |    

FEATURES...PERFORMANCE... ...AND TRAVEL READY!



The Go Go Elite Traveller® Plus offers easy disassembly and includes an ergonomic tiller for one-handed operation. Today's fast-paced world is no match for this stylish and portable scooter. With a 300 lb. weight capacity, a charger port in the tiller, one-hand feather-touch disassembly and maximum speeds up to 4.5 mph, it's easy to go where you want to go!

Additional Features

- Delta tiller with charger port to conveniently charge your scooter
- Convenient dual voltage charger can charge battery pack on-board or off-board
- Feather-touch disassembly permits simple frame separation with only one hand and easily disassembles into five super lightweight pieces
- Front frame-mounted seat post offers maximum stability
- Microprocessor-based controller offers optimal power management and added safety features
- Pride's exclusive black, non-scuffing tires
- Removable footboard basket for additional storage (4-wheel version only)

Colors (both color panel sets included)



Accessories*

- Cane clip attachment
- Cane/crutch holder (single, double)
- Cell phone holder
- Cup holder
- Forearm crutch holder
- 50" Lap belt
- Oxygen tank holder
- Rear basket
- Rearview mirror
- Saddlebag
- Safety flag
- Under seat storage (assembly required)
- Walker holder
- Weather cover
- Wishbone crutch holder
- XLR USB charger
- Extra 18AH battery box

Military Patches



* Only one rear accessory can be used at a time.



Two sets of easily changeable, red and blue colored panels.



Feather-touch disassembly for quick transport or storage.



Charger port located conveniently on the tiller.