

<b>Model</b>	Go Go® Endurance Li <sup>1</sup>
<b>Part Number</b>	S54LXLIT
<b>Weight Capacity</b>	325 lbs.
<b>Maximum Speed</b> <sup>2</sup>	Up to 4.2 mph
<b>Ground Clearance</b> <sup>4</sup>	3.25" at side deck
<b>Turning Radius</b> <sup>4</sup>	48.25"
<b>Overall Length</b> <sup>4</sup>	42"
<b>Overall Width</b> <sup>4</sup>	20.5"
<b>Seat-to-Deck Height</b> <sup>4</sup>	16.25" - 17.75"
<b>Seat-to-Ground Height</b> <sup>4</sup>	21.25"- 22.75"
<b>Front Tires</b> <sup>4</sup>	2" x 8" solid
<b>Rear Tires</b> <sup>4</sup>	2.5" x 8" solid
<b>Range Per Charge</b> <sup>2,3</sup>	Up to 7.2 miles (8AH)** Up to 14.4 miles (16AH MAXX) <sup>5</sup>
<b>Total Weight Without Batteries</b> <sup>4</sup>	93.5 lbs. w/o battery pack; with 8AH battery pack**: 101.5 lbs.; with 16AH MAXX battery pack <sup>5</sup> : 105.5 lbs.
<b>Heaviest Piece When Disassembled</b> <sup>4</sup>	40.75 lbs. (front section)
<b>Standard Seating</b>	Type: Mesh back Material: Black vinyl Dimensions: 17" x 17" Weight: 23 lbs.
<b>Motor Output</b>	24-volt DC motor, sealed transaxle
<b>Brakes</b>	Regenerative and electromechanical
<b>Battery Requirements</b> <sup>4</sup>	Size: 8 AH lithium ion pack** or 16 AH (MAXX) lithium ion pack <sup>5</sup> Weight: 8 lbs. (8AH); 12 lbs. (16AH MAXX)
<b>Battery Charger</b>	Off-board, 3.5 amp (lithium ion)
<b>Warranty</b>	Frame: Lifetime limited Electronics: 2-year limited Drive motors: 2-year limited Batteries: 2-year limited

## HCPCS Code: K0899

POWER MOBILITY DEVICE, NOT CODED  
BY DME PDAC OR DOES NOT MEET CRITERIA

- 1 Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- 2 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- 3 Results derived from theoretical calculations based on lithium battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 4 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- 5 Does not meet airline requirements.

**PRIDE**  
MOBILITY  
LIVE YOUR BEST®

(US) 800-800-8586 • (Canada) 888-570-1113  
pridemobility.com |      



\*\* Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Go Go® Endurance Li and its lithium ion battery are permitted.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

041624

**GO GO**<sup>®1</sup>  
ENDURANCE Li

**DISCOVER THE  
PORTABLE POWER  
OF LITHIUM ION**



FDA Class II Medical Device<sup>1</sup>

**LIVE YOUR BEST®**

pridemobility.com |      

# LITHIUM ION MEETS PORTABILITY



**LONGER  
LIFE**

## LITHIUM ION BATTERY

- Faster charging time
- Sustained power performance
- Lighter weight
- More useable energy
- Longer life
- Smart battery protection

Mesh seat or compact, foldable seat with armrests

Feather-touch disassembly



Delta tiller with charger port


High-intensity LED headlight

Interchangeable color panels for front and back

CTS Suspension for greater comfort



Charger port located conveniently in tiller.



Feather-touch disassembly for quick transport or storage.



Comfort-Trac front and rear suspension for added comfort.

With the **Go Go® Endurance Li**, traveling has never been easier! This 4-wheel scooter has everything you need for a safe and convenient trip. With the standard 8AH airline-compliant lithium-ion battery pack\*\*, the Go Go Endurance is perfect for air travel. Feather-touch disassembly provides the portability you need. When you arrive at your destination, enjoy a smooth, comfortable ride with Comfort-Trac Suspension.

## Additional Features

- MAXX per charge range up to 14.4 miles w/ 16AH battery pack<sup>5</sup>; 7.2 miles w/ 8AH battery pack\*\*
- Maximum speed up to 4.2 mph
- Delta tiller with ergonomic wraparound handles lets you operate the scooter with one hand and rest your wrists
- Convenient dual charging ports can charge battery pack on-board or off-board
- Includes two sets of easily changeable, colored panels
- Standard 8AH airline-compliant lithium-ion battery pack\*\* or upgrade to 16AH battery pack<sup>5</sup>
- Two seat options: 17"x17" mesh back or 20"x18" compact, foldable

\*\*Airline-compliant. Please see note on back page for details.

## Colors



Garnet Red (matte)

Euro Grey

## Accessories\*

- Cane clip attachment
- Cane/crutch holders
- Cell phone holder
- Cup holder
- Lap belt (50", 60", 70")
- Oxygen tank holder
- Rear basket
- Rear view mirror
- Saddlebag
- Under seat storage
- Walker holder
- Weather cover
- Wishbone crutch holder
- XLR USB charger
- Additional spare 8AH lithium-ion battery pack\*\*
- Additional spare 16AH MAXX lithium-ion battery pack<sup>5</sup>

## Military Patches



\*Only one rear accessory can be used at a time.