

PRIDE[®]
POWER
RECLINER
COLLECTION



FDA Class II
Medical Device¹

“

I LOVE MY PRIDE LIFT CHAIR!
I CAN FULLY RECLINE AND BE
COMPLETELY COMFORTABLE!
AND, THE FABRIC MATCHES
MY LIVING ROOM PERFECTLY!

– GEORGE R.



SCOTT MEUSER
CHAIRMAN & CEO
PRIDE MOBILITY®

With over three decades in the power lift recliner industry, we know that our business is about you, our customer. We're humbled when you choose to bring one of our power lift recliners into your home, knowing that it must fit your lifestyle and decor.

Pride® Power Lift Recliners is the culmination of your direct feedback. You wished for luxurious comfort, exceptional quality, state-of-the-art features and elegant designs. Over a 30-year period, we listened and refined, then listened some more. The result is the Pride line of power lift recliners designed for you, your home and your lifestyle.

We understand that bringing a power lift recliner into your home isn't just a purchase, but an investment – toward your comfort, home décor, and lifestyle. We deeply value that you, our customer, have enabled us to meet those vital goals through Pride Power Lift Recliners.




A handwritten signature in black ink that reads "Scott Meuser". The signature is fluid and cursive, with a long horizontal line extending from the end of the name.

RECLINE POSITION



UPWARD LIFT

FUNCTIONS AND FEATURES

-  1. The full recline position is ideal for napping or relaxing.
-  2. The upward lift offers convenient stand assist for getting in and out of the recliner.
-  3. Dual side pockets are ideal for storing remotes, books or magazines.




SIDE POCKETS

FOOTREST EXTENSION⁵



USB REMOTE



4. The footrest extension allows you to relax and stretch out completely.



5. The USB remote adjusts all functions to help you find the perfect position.



6. The backrest moves up and down independently to help you find the perfect position.



“

MY PRIDE RECLINER GIVES ME
THE SUPPORT I NEED
AND SO MUCH MORE!
I LOVE HOW EASY IT IS
TO ADJUST MY POSITION!

– EMILY R.


HIGH-RESILIENCY FIBER FILL | ■
for maximum comfort

FABRIC | ■
Soft, comfortable, luxurious fabrics are durable as they are stylish. Tested to commercial standards, guaranteed to last.



UNIQUE SPRING DESIGN | ■
Interlocking springs ensure maximum support while delivering individualized comfort.

STEEL FRAME | ■
Heavy-duty metal framework is forged to last.



STRONG WOOD FRAME
Premium furniture-grade hardwood is robust, durable and long-lasting.

MORTISE TENON JOINTS
Interlocking puzzle-like joints

HEAVY-DUTY LIFT ACTUATOR AND SCISSOR MECHANISM

Design also includes integrated battery backup and quiet and smooth lift system.

1. Pride® Power Recliner Collection products are FDA Class II Medical Devices designed to aid individuals with mobility impairments.



HERITAGE¹

FDA Class II Medical Device¹





There is a chill in the morning air as a new day begins. Before taking on the bustle of your life, embrace the inviting warmth and coziness of the Heritage collection by Pride® Power Lift Recliners. Stretch out completely with the footrest extension* and enjoy the steaming aroma from your morning coffee. Total Comfort Seating ensures outstanding relaxation and comfort, as you feel the softness of luxurious Cloud 9 fabrics enveloping your body.

When you can't put the day off any longer, the Heritage Collection lifts you quickly and smoothly, helping you to get up and moving. With a quiet and smooth lift system that offers durable and long-lasting performance, the Heritage collection fulfills your every need. For the recliner that perfectly complements every aspect of your life, choose the Heritage collection by Pride.



STANDARD FABRIC: CLOUD 9



BLACK CHERRY PACIFIC STONE WALNUT

*Footrest extension is not available on the LC-358S, LC-358XL & LC-358XXL.





LC-105

When it comes to the essentials of everyday life, never sacrifice your comfort. Experience a simply elegant solution with the Essential LC-105 Power Lift Recliner. Fill the empty spot in your house or apartment with the chair that complements your décor perfectly. With beautiful micro-suede fabrics in soft, dark colors, the Essential Lift Chair from Pride® Power Lift Recliners is the accent piece you've been waiting for!

The Essential LC-105 has an easy-to-use hand remote, so finding a comfortable position is as simple as pressing a button! Sinuous Spring Seating gives you endless support. Rolled, pleated armrests add a touch of elegance. With dual pockets that are easily accessible, the essentials of your life are always at your fingertips. Discover the elegant simplicity of the Essential LC-105 Recliner.



STANDARD FABRIC: MICRO-SUEDE




OPTIONAL FABRIC:
LEXIS URETHANE®



COCOA



SANDAL



NEW
CHESTNUT

ESSENTIAL¹



LC-250

Although the demands of life are relentless, take solace in the comfort of a Pride® Power Lift Recliner. The Essential LC-250 Recliner gives you the boost you need to take on the day. Total Comfort Seating provides complete support for ultimate relaxation. The two-button hand remote is simple to use. Adjust your chair effortlessly to find your perfect position every time.

With luxurious Cloud 9 fabrics that are soft to the touch, you'll find the perfect match to suit your home décor. Enjoy the plentiful storage of dual pockets that are easy to reach. Experience quiet and smooth lift performance, day after day, year after year. It's time to invest in your comfort with the Essential LC-250 Recliner.



STANDARD FABRIC: CLOUD 9



PACIFIC



STONE



WALNUT



CHESTNUT

OPTIONAL FABRIC:
LEXIS STA-KLEEN®

ESSENTIAL¹

OPTIONAL FABRICS

ULTRAFABRICS® ULTRALEATHER®



BUFF FUDGE GARNET PECAN CHARCOAL

Ultrafabrics® Ultraleather® on Pride® Power Lift Recliners offer a richly detailed grain pattern expressive of premium grade leather. Ultrafabrics polyurethane products are manufactured utilizing Tukumi™ technology, a 4-layer laboratory process to create a luxurious, durable and comfortable polyurethane fabric.

FEATURES AND BENEFITS

- Rich feel and natural grain pattern expressive of soft luxury-grade leather
- Clinically proven more resistant to temperature changes than vinyl and genuine leather for comfort
- Antimicrobial surface cleans effortlessly with soap and water

CRYPTON® SUPER FABRIC




ESPRESSO COOL GREY RED SAND

Crypton® Super Fabrics from Pride® Power Lift Recliners are not only luxurious, durable and contemporary, but also resist spills, stains, odors and bacteria. Nothing gets past Crypton's unique patented integrated fiber process, protecting your upholstery inside and out.

FEATURES AND BENEFITS

- 100% water resistant – spills just wipe away
- Only disinfectable fabrics in the U.S. approved by the EPA
- Permanent technology that will not come off the fabrics
- 5-year warranty on moisture barrier

LEXIS STA-KLEEN® PERFORMANCE FABRIC



CHESTNUT MUSHROOM BLACK BURGUNDY

Lexis Sta-Kleen® Performance Fabric is a revolution in faux leather upholstery — the first stain resistant urethane fabric. Sta-Kleen offers permanent protection, created with a proprietary bonding process inherent in the manufacturing process of the product. It is not a topical application; it will last the lifetime of the upholstery.

FEATURES AND BENEFITS

- Look and feel of real leather, but with care-free durability
- Exceptional strength, abrasion resistance, and UV stability
- "Stain shield" allows spills to wipe off with a cloth, and it's anti-bacterial

QUICK REFERENCE FABRIC GUIDE

MODEL	CLOUD 9	MICRO-SUEDE	LEXIS STA-KLEEN	CRYPTON ARIA	ULTRALEATHER	LEXIS URETHANE
LC-358PW	Standard	N/A	Optional	Optional	Optional	N/A
LC-358S	Standard	N/A	Optional	Optional	Optional	N/A
LC-358M	Standard	N/A	Optional	Optional	Optional	N/A
LC-358L	Standard	N/A	Optional	Optional	Optional	N/A
LC-358XL	Standard	N/A	Optional	Optional	Optional	N/A
LC-358XXL	Standard (only Stone & Walnut)	N/A	Optional (Chestnut only)	N/A	N/A	N/A
LC-105	N/A	Standard	N/A	N/A	N/A	Optional
LC-250	Standard	N/A	Optional (Chestnut only)	N/A	N/A	N/A

ACCESSORIES



HEAD & LUMBAR PILLOW

(SOLD SEPARATELY)⁵

Adds additional support and comfort (in the same fabric the chair was ordered in) in the head or lumbar area.



BOLSTER PILLOW

(AVAILABLE ON MOST PRIDE® POWER LIFT RECLINERS)⁵

Fabrics to match your chair. Round, firm pillow with zipper to adjust the fill. Adds additional customized support where needed.



NON-SKID LEG LEVELERS

Non-skid leg levelers keep your Power Lift Recliner from sliding on a vast array of flooring.



HEAT AND MASSAGE⁴

Heat and massage offers two massage units and one heat pad.



CLEANING KIT

This kit can be used with any Power Lift Recliner fabric (excluding vinyl/leather) in order to remove almost any stain or odor.



LITHIUM BATTERY BACKUP

(NOT AVAILABLE ON LC-105 & LC-250)

Next time the electricity goes out, be prepared with the new optional rechargeable lithium backup battery.

SPECIFICATIONS



LC-105



LC-250



LC-358

COLLECTION	ESSENTIAL ¹	ESSENTIAL ¹
MODEL	LC-105	LC-250
POSITION	3-Position 5'3" & Below	3-Position 5'4" - 5'9"
SIZE	Small	Medium
WEIGHT CAPACITY	325 lbs.	375 lbs.
SEAT-TO-FLOOR	18"	20"
SEAT DEPTH	19"	20"
SEAT WIDTH	20"	22"
TOP OF BACK TO SEAT	24"	26"
FULLY PADDED CHAISE	Yes	Yes
BACK STYLE	Sewn Pillow Back	Blown Fiber Split Back
ADDITIONAL OPTIONS	No	No
FOOTREST EXTENSION	N/A	N/A
SEATING TYPE	Sinuous Spring	Total Comfort
DISTANCE FROM WALL	17"	20"
OVERALL WIDTH	33"	34"
CHAIR WEIGHT	96 lbs.	118 lbs.
WARRANTY	1-Yr Limited	2-Yr. Limited

1. Pride® Power Recliner Collection products are FDA Class II Medical Devices designed to aid individuals with mobility impairments.
2. Available on Heritage Collection except LC-358XL & XXL
3. Available on LC-358 S, PW, M & L
4. USB hand control not available with heat & massage options.
5. Not available on all models.

* Not all units achieve all positions.
Please refer to specific product.

COLLECTION	HERITAGE ¹	HERITAGE ¹	HERITAGE ¹	HERITAGE ¹	HERITAGE ¹	HERITAGE ¹
MODEL	LC-358PW	LC-358S	LC-358M	LC-358L	LC-358XL	LC-358XXL
POSITION	3-Position 5'3" & Below	3-Position 5'3" & Below	3-Position 5'4" - 5'9"	3-Position 5'9" - 6'	3-Position 5'8" - 6'1"	2-Position 5'8" - 6'1"
SIZE	Petite Wide	Small	Medium	Large	Heavy Duty	Very Heavy Duty
WEIGHT CAPACITY	375 lbs.	375 lbs.	375 lbs.	375 lbs.	500 lbs.	600 lbs.
SEAT-TO-FLOOR	19"	19"	20"	21"	20"	20"
SEAT DEPTH	19"	19"	20"	22"	22"	22"
SEAT WIDTH	24"	18"	20"	22"	26"	30"
TOP OF BACK TO SEAT	23.5"	24"	24.5"	25.5"	29"	29"
FULLY PADDED CHAISE	Yes	Yes	Yes	Yes	Yes (removable)	Yes (removable)
BACK STYLE	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back
ADDITIONAL OPTIONS	Heat & Massage ⁴ Lithium Battery Backup	Heat & Massage ⁴ Lithium Battery Backup	Heat & Massage ⁴ Lithium Battery Backup	Heat & Massage ⁴ Lithium Battery Backup	Heat & Massage ⁴ Lithium Battery Backup	Lithium Battery Backup
FOOTREST EXTENSION	Standard	N/A	Standard	Standard	N/A	N/A
SEATING TYPE	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Sinuous Spring	Sinuous Spring
DISTANCE FROM WALL	18"	18"	17"	19"	19"	14"
OVERALL WIDTH	40"	35"	35"	37"	39.5"	43"
CHAIR WEIGHT	130 lbs.	118 lbs.	117 lbs.	120 lbs.	159 lbs.	168 lbs.
WARRANTY	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; We reserve the right to alter specifications without prior notice.



Look for the Quick Ship Indicator for the shortest lead times available.

PRIDE[®]
POWER
RECLINER
COLLECTION

(US) 800.800.8586 | (CAN) 888.570.1113 | pridemobility.com |

