



PRIDE[®]
M O B I L I T Y
L I V E Y O U R B E S T[®]

TRAVEL/SCOOTER COLLECTION

4-5
PX4

6-7
Pursuit 2

8-9
Zero Turn 8 / Zero Turn 10

10-11
iRIDE® 2 / i-Go™

12-13
Go Go Elite Traveller® / Endurance Li

14-15
Go Go® Sport / LX with CTS

16-17
Go Go® Ultra X / ES 2

18-19
Revo® 2.0 / Victory® LX Sport

20-21
Victory® 10 / 10.2

22-23
Victory® 9 / Maxima

24
Accessories



“

MY PRIDE SCOOTER HAS
BEEN A LIFESAVER
FOR ME. I COULD NOT
DO LIFE WITHOUT IT.

— Etta G.,
White Plains, NY

PX4

The PX4 features an impressive 500-lb. weight capacity, a top speed of 9.3 mph, water-resistant user interface, enhanced lighting package, and a remarkable 23.4 mile battery range.

The ignition system also offers an extra layer of security with its innovative NFC technology. Simply tap or place the NFC key fob near the receiver to start the scooter.



Up To
23.4
Miles per Charge

82.3"
Turning Radius

Up To
9.2mph
Maximum Speed

500lbs.
Weight Capacity

FDA Class II Medical Device¹

Forest Camo (Limited Edition)

Satin Aluminum

Dark Violet

Red

Peacock Blue

Sunflower



MORE ROOM TO ROAM



PURSUIT²

Pursuit 2

Built to thrill and designed to go the distance, the Pursuit 2 never backs down from an off-road detour or long-distance excursion. With a fast top cruising speed of 9.3 mph and the ability to make it 40 miles on a single charge, this agile unit offers swift, efficient transport and endurance you can count on.

You were born to roam. The Pursuit 2 was built to ride.



Up To
40
Miles Per Charge*

69.3"
Turning Radius

Up To
9.3mph
Maximum Speed

400lbs.
Weight Capacity

FDA Class II Medical Device¹

* Requires dual lithium battery upgrade

Grey (matte)

Orange

Black



INSERT YOURSELF: ANYWHERE



ZERO TURN8

Zero Turn 8

Your mobility has turned a corner with the Zero Turn 8, an innovative scooter equipped with iTurn Technology™, enabling you to maneuver this scooter through small spaces and around tight corners! Plus, enjoy full lighting, plenty of storage and a USB charger.



Up To
13.1
Miles Per Charge at 200 lbs.

37.25"
Turning Radius

Up To
5.9mph
Maximum Speed

325lbs.
Weight Capacity

FDA Class II Medical Device*

Sapphire Blue

Pewter

Pearl White

Black Onyx

Green Moonstone

Fire Opal



ZERO TURN10

Up To
21.2
Miles Per Charge at 200 lbs.

43"
Turning Radius

Up To
7.3mph
Maximum Speed

400lbs.
Weight Capacity



Zero Turn 10

The Zero Turn 10 can traverse the most beautiful, out-of-reach places and move seamlessly through doors and hallways in your home. This scooter combines powerful dual motors and two-wheel drive with revolutionary turning technology. Let the Zero Turn 10 carry you to places unknown.

FDA Class II Medical Device¹

Candy Apple Red

Ocean Blue (Matte)

Pearl White

Black¹¹ (Matte)



iRIDE²

iRide[®] 2

Journey far and wide with the iRIDE[®] 2! This ultra-lightweight scooter is super portable and perfect for traveling. Equipped with both a manual handbrake and supplemental foot brake, the iRIDE[®]2 delivers a safe and comfortable ride. A folding back has been added for increased comfort, along with an upgraded LCD display for increased functionality.



Up To
12
Miles Per Charge at 200 lbs.

33.5"
Turning Radius

Up To
5.8mph
Maximum Speed

250lbs.
Weight Capacity

Non-Medical Device

Raspberry

Mango

Arctic Ice

Blackberry





Up To
9.3
Miles Per Charge at 220 lbs.

45.9"
Turning Radius

Up To
4mph
Maximum Speed


275lbs.
Weight Capacity



Non-Medical Device

i-Go™

The i-Go™ Folding Scooter supports you when you refuse to compromise. So, when you're dreaming of a cruise vacation, know that i-Go™ comes with the all-day power, speed and versatility to take you there and back again. When disassembled, the heaviest piece weighs less than 32 lbs. and has a height of less than 12" for easy storage.



FLIP-UP ARMRESTS
(OPTIONAL ACCESSORY)



PHONE HOLDER
(OPTIONAL ACCESSORY)



REMOVABLE STORAGE
(STANDARD)



- Iceberg Blue (matte)
- Sugar Plum (matte)
- White (matte)
- Black (matte)

- ACCESSORIES:**
- Armrests
 - Phone holder
 - Cup holder
 - Rearview mirror
 - Spare battery
 - Travel kit

GO GO
ELITE TRAVELLER®

Go Go Elite Traveller®

Every day is a new adventure with the Go Go Elite Traveller®, an exceptionally stable and maneuverable travel scooter. With easy disassembly and a convenient storage basket, the Go Go Elite Traveller keeps up with your active life!



Up To
9.75
Miles Per Charge

44"
Turning Radius

Up To
4mph
Maximum Speed

300lbs.
Weight Capacity

FDA Class II Medical Device¹

Also available in 3-wheel. Please see back two pages for specifications.

Red

Blue

Interchangeable Shrouds



Up To
14.4
Miles Per Charge
(16AH* MAXX Battery Pack)

48.25"
Turning Radius

Up To
4.2mph
Maximum Speed

325lbs.
Weight Capacity

GO GO
ENDURANCE Li



Go Go® Endurance Li

With the Go Go® Endurance Li, traveling has never been easier! This 4-wheel scooter has everything you need for a safe and convenient trip. With the standard 8AH airline-compliant lithium battery pack, the Go Go Endurance is perfect for air travel. Feather-touch disassembly provides the portability you need.

*16 AH battery pack does not meet airline requirements.

FDA Class II Medical Device¹

Garnet Red (matte)

Euro Grey

Interchangeable Shrouds



GOGO
SPORT

Go Go® Sport

It's time to make your mark on this world! The Go Go® Sport offers simple disassembly and an ergonomic delta tiller, perfect for one-handed operation. Enjoy the versatility of this sporty scooter, as well as bright LED lighting.



Up To
12
Miles Per Charge

50.5"
Turning Radius

Up To
4.7 mph
Maximum Speed

325 lbs.
Weight Capacity

FDA Class II Medical Device¹

Also available in 3-wheel. Please see back two pages for specifications.

Red

Blue



GO GO[®]
LX WITH CTS SUSPENSION

Up To
6
Miles Per Charge

33.5"
Turning Radius

Up To
4.4mph
Maximum Speed

300lbs.
Weight Capacity



Go Go[®] LX with CTS

There's smooth, and then there's the Go Go[®] LX with CTS! This scooter has revolutionary suspension for the smoothest, most comfortable ride. From the park to your backyard, this feature-rich scooter offers best-in-class ride quality.

FDA Class II Medical Device¹



Red

Blue

Interchangeable Shrouds

Orion Yellow

Onyx Black

Pearl Pink

Sunburst Orange

Optional



GO GO[®]
ULTRA X

Go Go[®] Ultra X

Enjoy worry-free travel with the Go Go[®] Ultra X, an ultra-compact scooter that is simple to disassemble. Ride in comfort and style with the foldable swivel seat and keep your personal items handy with the front basket.



Up To
8.9
Miles Per Charge at 200 lbs.

52.25"
Turning Radius

Up To
4.1mph
Maximum Speed

300lbs.
Weight Capacity

FDA Class II Medical Device¹

Also available in 3-wheel. Please see back two pages for specifications.

Black

PRIDE



GO•GO®
es2

Up To
6.24
Miles Per Charge

31"
Turning Radius

Up To
4mph
Maximum Speed

250lbs.
Weight Capacity



Go Go® ES 2

Get ready to go with the Go Go® ES 2 Scooter! Compact and lightweight, this scooter has two speed settings for ultimate versatility, whether you are indoors or outdoors. It's perfect for errands and quick trips! Plus, a front basket provides plenty of storage space.

FDA Class II Medical Device¹

Black



REVO® 2.0

Revo® 2.0

Robustly rugged meets brilliant convenience. With the Revo® 2.0, there is nothing you can't handle! This tough and dependable scooter has powerful suspension for a smooth ride, yet easily disassembles for full portability.



Up To
17.8
Miles Per Charge at 200 lbs.

52"
Turning Radius

Up To
5.2mph
Maximum Speed

400lbs.
Weight Capacity

FDA Class II Medical Device¹

Also available in 3-wheel. Please see back two pages for specifications.

Grey Street

True Blue



VICTORY® LX^{SPORT}

Up To
19.7
Miles Per Charge at 200 lbs.

61.25"
Turning Radius

Up To
8.6mph
Maximum Speed

400lbs.
Weight Capacity



FDA Class II Medical Device¹

Victory® LX Sport

You're active and always on the go. You need a scooter that is up to the task! Appeal to your sporty side with the Victory® LX Sport. Enjoy premium suspension, bright lighting and other sophisticated features, all in one sporty package.

Candy Apple Red

True Blue



VICTORY 10

Victory® 10

Versatility knows no bounds with the Victory® 10. Enjoy simple disassembly with one hand and drive up to 16 miles on a single charge. You won't miss a beat with this scooter, which includes a bright LED headlight. Anything is possible with the Victory 10!



Up To
15.5
Miles Per Charge

54"
Turning Radius

Up To
5.3mph
Maximum Speed

400lbs.
Weight Capacity

FDA Class II Medical Device¹

Also available in 3-wheel. Please see back two pages for specifications.

Candy Apple Red

Viper Blue



VICTORY 10.2

Up To
13
Miles Per Charge

45"
Turning Radius

Up To
5.2mph
Maximum Speed

400lbs.
Weight Capacity



Victory® 10.2

Embrace your potential! The Victory® 10.2 is fully loaded with style and advanced features, including frosted LED headlights and an integrated storage space with cup holder. Drive at speeds up to 5.2 mph and see what life has in store for you!

ONLY available in 3-wheel. Please see back two pages for specs.

FDA Class II Medical Device¹

Candy Apple Red

Ocean Blue (Matte)



VICTORY® 9

Victory® 9

Sophisticated and stylish, the Victory® 9 is everything you could want in a scooter. The delta tiller lets you operate the scooter with a single hand. Enjoy easy disassembly and plenty of lighting with the LED headlight.



Up To
13
Miles Per Charge

51.5"
Turning Radius

Up To
5.3mph
Maximum Speed

300lbs.
Weight Capacity

FDA Class II Medical Device¹

Also available in 3-wheel. Please see back two pages for specs.

Candy Apple Red



MAXIMA

Up To
19.5
Miles Per Charge

53"
Turning Radius

Up To
5.8mph
Maximum Speed

500lbs.
Weight Capacity



Maxima

Experience the comfort of a new high-back, depth adjustable limited recline seat that swivels for easier access.

The Maxima scooter comes standard with front and side bumpers, a USB charger and full LED lighting!

Also available in 3-wheel. Please see back two pages for specs.


FDA Class II Medical Device¹

Garnet Red (Matte)

Ocean Blue (Matte)



ACCESSORIES



1. Cane Clip Attachment | 2. Cup Holder | 3. Under Seat Storage | 4. Single OR Double Cane Holder | 5. XLR USB CHARGER
6. Rear Basket | 7. Oxygen Holder | 8. Ram® X-Grip Cell Phone Holder | 9. Wishbone Crutch Holder | 10. Lap Belt | 11. Saddlebag
12. Walker Holder | 13. Rearview Mirror | 14. Weather Cover | 15. Military Patches | 16. Safety Flag

“







MY PRIDE SCOOTER IS
SIMPLE TO OPERATE
AND A BREEZE TO TRANSPORT.
I TAKE MY SCOOTER EVERYWHERE
WITH ME. PRIDE SCOOTERS HAVE
GIVEN ME MY LIFE BACK!

— Sandy M.,
Homer Glen, IL

MODEL NAME / NUMBER	PX4 ¹ SC134	PURSUIT 2' SC7132	ZERO TURN 8' JAZZYZTS	ZERO TURN 10' S710ZT	iRIDE® 2 S25	i-GO™ SC20	GO GO ELITE TRAVELLER® ¹ SC40E (3W), SC44E (4W)	GO GO® ENDURANCE Li ¹ S54XLIT	GO GO® SPORT ¹ S73 (3W), S74 (4W)
WEIGHT CAPACITY	500 lbs.	400 lbs.	325 lbs.	400 lbs.	250 lbs.	275 lbs.	300 lbs.	325 lbs.	325 lbs.
MAXIMUM SPEED ²	Up to 9.2 mph	Up to 9.3 mph	Up to 5.9 mph	Up to 7.3 mph	Up to 5.8 mph	Up to 4 mph	Up to 4 mph	Up to 4.2 mph	Up to 4.7 mph
GROUND CLEARANCE ^{3&8,9}	3.5" at motor frame, 4.4" at mid-frame	3.78" at transaxle, 4.6" at mid-frame	1.3" at motor, 1.5" at anti-tip, 3.5" (side deck)	2" at motor 3" at mid-frame	1.3" at mid-frame	4.2" at mid-frame	1.38" at motor, 0.75" at anti-tip, 2.25" at mid frame	1.5" at motor 3.23" at deck	1.5" at motor, 0.75" at anti-tip
TURNING RADIUS ³	82.3"	69.3"	37.25"	43"	33.5"	45.9"	3W: 33" 4W: 44"	48.25"	3W: 36" 4W: 50.5"
OVERALL LENGTH / WIDTH ^{3,7}	56" / 27-27.8"	57.5" / 29.7"	45.75" / 20.75"	48" / 24.25"	38" / 19.25"	44" / 20.3"	3W: 37" 4-W: 39.5" / 19.5"	42" / 20.5"	3W: 42" 4W: 43" / 21.25"
SEAT-TO-GROUND HEIGHT	24.17-25.29"	28.74"	22" - 23.5"	23.25" - 24.75"	19" - 20"	22.25"	18.25"-20.25"	21.25"-22.75"	20.5" - 22"
SEAT-TO-DECK HEIGHT	16.54-18.39"	17.6-20.9"	17" - 18.5"	17.75"-19.25"	11.5" - 12.5"	14.6"	15.75"-17.25"	16.25"-17.75"	17" - 18"
FRONT WHEELS (TIRES)	13" pneumatic	13" pneumatic	2"x7" solid	9" solid	6.5" solid	8" hub motor	3W: 2" x 8" solid 4W: 2" x 7" solid	2" x 8" solid	3W: 3" x 8" solid 4W: 3" x 9" solid
REAR WHEELS (TIRES)	13" pneumatic	13" pneumatic	2.5"x9" solid	10.75" solid	6.75" solid	9" perforated	2.5" x 8" solid	2" x 8" solid	3" x 9" solid
SUSPENSION	Front & Rear	Front & Rear	Front and rear (CTS)	Front and rear (CTS)	N/A	N/A	N/A	Front and rear (CTS)	N/A
RANGE PER CHARGE ⁵	Up to 23.4 miles (220 lbs.) Up to 14 miles at (500 lbs.)	Single Battery: Up to 17.4 miles (220 lbs.) Up to 13.8 miles (400 lbs.) Dual Battery: Up to 34.8 miles (220 lbs.) Up to 30.9 miles (350 lbs.)	Up to 11 miles (at 325 lbs.) Up to 13.1 miles (at 200 lbs.)	Up to 15.9 miles (at 400 lbs.) Up to 21.2 miles (at 200 lbs.)	Up to 10 miles at 250 lbs. Up to 12 miles at 200 lbs.	Up to 7 miles at 275 lbs. Up to 9.3 miles at 220 lbs.	3W: Up to 7.25 miles (12 AH) 4W: Up to 6 miles (12 AH) 3W: Up to 10.5 miles (18 AH) 4W: Up to 9.75 miles (18 AH)	8 AH: Up to 7.2 miles 16 AH: Up to 14.4 miles at 200 lbs.	3W: Up to 12.75 miles 4W: Up to 12 miles
HEAVIEST PIECE ³	201.2 lbs. (base)	163 lbs. (base)	48 lbs. (rear section)	75 lbs. (front section)	32 lbs. (base only)	32.8 lbs. (rear base with seat back)	3W: 29.5 lbs. (front section) 4W: 34 lbs. (front section)	40.75 lbs. (front section)	3W: 35.75 lbs. (front section) 4W: 44 lbs. (front section)
STANDARD SEATING	Limited Recline High Back (18" x 18")	Limited Recline High Back with Sliders (18" x 18")	Type: Compact, foldable Weight: 27 lbs. Material: Black vinyl Dimensions: Width: 17", Depth: 17", Back height: 14"	Type: Foldable, deluxe contour Weight: 34 lbs. or Foldable, compact Weight: 30 lbs. Material: Black vinyl Dimensions: Width: 18", Depth: 17"	Type: Compact, removable with folding back Weight: 10 lbs. Material: Molded plastic Dimensions: 14.75" x 13.5"	Type: Compact, folding (back is removable) Weight: 10 lbs. Material: Molded plastic Dimensions: 16" x 16"	Type: Compact, foldable Weight: 23 lbs. Material: Black vinyl Dimensions: Width: 17", Depth: 17"	Type: Compact, foldable Weight: 23 lbs. Material: Black vinyl Dimensions: 17"x17", 20"x18"	Type: Compact, foldable Weight: 23 lbs. Material: Black vinyl Dimensions: Width: 17", Depth: 17"
BATTERIES ^{4,13}	(2) 12V 75Ah (SLA/AGM)	(1) or (2) 24V, 50Ah Lithium-ion	(2) 22 AH	(2) 40 AH	(1) 36 volt, 5.2 AH lithium ion ¹³	(1) 25.2V 9.6 Ah Rechargeable Lithium-Ion	(2) 12 AH (standard) or (2) 18 AH (optional)	8 AH ¹³ lithium ion battery pack or 16 AH ¹² lithium ion battery pack (optional)	(2) 18 AH
BATTERY WEIGHT ^{5,6}	48.8 lbs. each	23.8 lbs. each	32 lbs. (battery box), 13.5 lbs. each	30.25 lbs. each	4.5 lbs.	4.8 lbs.	9 lbs. each (12 AH), 20 lbs. per battery pack 14 lbs. each (18 AH), 29.5 lbs. per battery pack	8 lbs. (8 AH); 12 lbs. (16 AH)	14 lbs. each 29.5 lbs. per battery pack
BATTERY CHARGER	Off-board, 8A	Off-board, 8A	Off-board, 3.5 amp	Off-board, 3.5 amp	36 volt lithium, off-board	Off-board, 3.5 amp (lithium-ion)	Off-board, 2 amp	Off-board, 3.5 amp (lithium ion) ¹³	Off-board, 2 amp
HCPCS CODES	K0899*	K0899*	K0899*	K0899*	N/M (non-medical)	N/M (non-medical)	K0800 (Group 1, Standard)	K0899*	K0801 (Group 1, Heavy Duty)
WARRANTY	Frame: Lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Batteries: 6-month limited	Frame: Lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Batteries: 6-month limited	Frame: Lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Batteries: 6-month limited	Frame: Lifetime limited Drive motors: 2-year limited Electronics: 2-year limited Batteries: 6-month limited	Frame: 3-year limited Electronics: 2-year limited Drivetrain, motor, battery and charger: 1-year limited	Frame: 3 years Electronics: 2 years (does not include charger); 1 year on motor, battery and charger	Frame: Lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Batteries: 6-month limited	Frame: Lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Batteries: 2-year limited	Frame: Lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Batteries: 6-month limited

PRIDE®

MOBILITY
LIVE YOUR BEST®

(US) 800-800-8586 | (CAN) 888-570-1113 | pridemobility.com |      

1. Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
2. Varies with user weight, terrain type, battery amp- hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.
3. Varies with user weight, terrain type, battery amp- hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.
4. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%.
5. Deep cycle (AGM or gel cell type recommended).
6. Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
7. Battery weight may vary based on manufacturer.
8. Width of base (not including seat).
9. Does not account for anti-tip clearance.
10. At midframe.
11. Product is available exclusively to Authorized Pride Brick & Mortar Providers. No internet sales are permitted. Available in store only. Other conditions apply.
12. Does not meet airline requirements.
13. Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your item and its lithium ion battery are permitted.

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.

N/M: This product is not a medical device and is not intended to assist, treat, diagnose, or alleviate any medical condition or disability.

*K0899 – power mobility device, not coded by DME PDAC.

**Not all accessories are available for all products. Please speak with your Provider for details and availability.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

