

SOPHISTICATED AND STYLISH

VICTORY^{® 1}

Versatility knows no bounds with Victory[®] Scooters.¹ Enjoy simple disassembly with one hand and advanced features designed for convenience and peace of mind. You can feel confident that the Victory's long-lasting batteries will keep you going. Anything is possible with a Victory[®] Scooter!

All models feature...

- Feather-touch disassembly
- Standard front basket
- Delta tiller with port to conveniently charge your scooter



FDA Class II Medical Device¹



SCAN THE QR CODE
TO LEARN MORE ABOUT
VICTORY[®] SCOOTERS!



FOLLOW US ON SOCIAL MEDIA

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

VICTORY® 10.2¹



FDA
Class II
Medical
Device¹

SPEED²
Up to
5.2 mph
—
**WEIGHT
CAPACITY**
400 lbs.
—
COLORS

FEATURES

- Available in a 3-wheel model only
- Frosted LED headlights and low-mounted LED pathway light brightens your route
- 18"x17" seat with swivel and sliders
- Integrated storage/cup holder in tiller
- Maximum range up to 13 miles (40AH batteries)^{3,4}

WARRANTY

Lifetime limited on frame; 2-year limited on electronics
2-year limited on drivetrain; 6-month limited on batteries

VICTORY® 10¹



FDA
Class II
Medical
Device¹

SPEED²
Up to
5.3 mph
—
**WEIGHT
CAPACITY**
400 lbs.
—
COLORS

FEATURES

- Available in 3 and 4-wheel models
- LED headlight
- 18"x17" black vinyl seat with 8-point swivel and sliders
- Ambient voltmeter
- Maximum range up to 15.5 miles with 3 and 4-wheel models³

WARRANTY

Lifetime limited on frame; 2-year limited on electronics
2-year limited on drivetrain; 6-month limited on batteries

VICTORY® 9¹



FDA
Class II
Medical
Device¹

SPEED²
Up to
5.3 mph
—
**WEIGHT
CAPACITY**
300 lbs.
—
COLORS

FEATURES

- Available in 3 and 4-wheel models
- LED headlight
- 18"x17" black vinyl seat with 8-point swivel and sliders
- Ambient voltmeter
- Maximum range up to 13 miles³

WARRANTY

Lifetime limited on frame; 2-year limited on electronics
2-year limited on drivetrain; 6-month limited on batteries

K0899 - POWER MOBILITY DEVICE,
NOT CODED BY DME PDAC OR DOES
NOT MEET CRITERIA

VICTORY® LX SPORT¹



FDA
Class II
Medical
Device¹

SPEED²
Up to
8.6 mph⁵
—
**WEIGHT
CAPACITY**
400 lbs.
—
COLORS

FEATURES

- Available in a 4-wheel model only
- CTS Suspension
- Full LED lighting package includes curb light, rear lights and directional signals
- Puddle light under tiller
- 18"x18" high back seat with sliders
- Maximum range up to 19.7 miles with 200 lb. user weight and up to 13 miles with 400 lb. user weight³

WARRANTY

Lifetime limited on frame; 2-year limited on electronics
2-year limited on drivetrain; 6-month limited on batteries

1) Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
2) Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This can be subject to a variance of (+ or -) 10%.
3) Range varies based on user weight, terrain and battery type.
4) Maximum range up to 11 miles with U-1 batteries.
5) Maximum speed is reduced to 8 mph with a user weight over 200 lbs.