



**PRIDE**<sup>®</sup>  
MOBILITY  
LIVE YOUR BEST<sup>®</sup>

SIGNATURE COLLECTION



**4-7**

Go Go® Endurance Li

---

**8-11**

Zero Turn 10

---

**12-15**

Zero Turn 8

---

**16-19**

Jazzy Air® 2

---

**20-23**

Jazzy® Carbon

---

**24-27**

Go Chair®

---

**28-29**

Accessories

---

**30-31**

Specifications



**SCOTT MEUSER**

CHAIRMAN & CEO  
PRIDE MOBILITY®

With 37 years of experience under our belts, Pride realizes our business is about you. Our products are built with your input in mind and are designed to meet the needs of people in the real world. We offer a mobility product for every person and lifestyle.

The success of Pride Mobility Products is driven by our strong passion to make a difference in people's lives by providing innovative, market-leading products and friendly service. This begins with the dedication and loyalty of over 1,000 Pride employees: Their hard work, coupled with their skills, knowledge and expertise, supports our mission to be the friendliest company in the industry.

Together, with you, we want to continue our operational excellence by listening to your feedback and bringing it to life in our products and services. Thank you for allowing us into your homes and lives, and for being part of the Pride Family.

A handwritten signature in black ink that reads "Scott Meuser". The signature is fluid and cursive, with a long horizontal line extending to the right.

# THE LITHIUM ION ADVANTAGE

- SUSTAINED POWER PERFORMANCE
- MORE USEABLE ENERGY
- SMART BATTERY PROTECTION



FDA Class II Medical Device\*



16 AH MAXX batter box  
NOT AIRLINE COMPLIANT

## PORTABLE POWER. LITHIUM ION.



**LONGER**  
LIFE

**140%**  
**FASTER**  
CHARGING TIME

**60%**  
**LIGHTER**  
WEIGHT

**300%**  
**LONGER**  
BATTERY LIFE CYCLE

**SUSTAINED**  
**POWER PERFORMANCE**

**GO GO®**  
**ENDURANCE Li**

With the Go Go® Endurance Li, traveling has never been easier! This 4-wheel scooter has everything you need for a safe and convenient trip. With the standard 8AH airline-compliant lithium battery pack, the Go Go Endurance is perfect for air travel. Feather-touch disassembly provides the portability you need.

# GOGO<sup>®</sup> ENDURANCE Li

Garnet Red  
(matte)

Euro Grey

Compact, foldable  
seat with armrests

Feather-touch  
disassembly

CTS Suspension for greater comfort

Delta tiller with  
charger port

High-intensity  
LED headlight

Interchangeable  
color panels for  
front and back

FDA Class II Medical Device \*

**325lbs.**

Weight Capacity

Up To

**4.25mph**

Maximum Speed


**48.25"**

Turning Radius

Up To

**14.4**

Miles Per Charge with  
16AH MAXX Battery Pack







# REIMAGINE INDEPENDENCE

## **ZERO TURN10**

Imagine your perfect scooter. A sleek, contemporary design that can traverse the most beautiful, out-of-reach places. Imagine that same scooter in your home, moving seamlessly through doors and down hallways. Your imagination becomes reality with the **Zero Turn 10, (ZT10)** a scooter that combines powerful dual motors and two-wheel drive with revolutionary turning technology. Let the Zero Turn 10 carry you to places unknown.

# ZERO TURN10

- Black\*\*\*  
(matte)
- Candy Apple  
Red
- Ocean Blue  
(matte)
- Pearl White

Foldable, deluxe contour seat with swivel and sliders, height width, depth and angle-adjustable armrest

LED lights with sequential turn signals

New, improved ambient voltmeter

Puddle light

USB charger

Cupholder/storage

Tiller adjustment lever

Standard large basket

LED headlights with sequential turn signals

Comfort-Trac Suspension (CTS)

FDA Class II Medical Device \*



Our patent-pending, intelligent-turning *iTurn Technology™* ensures effortless navigation in small spaces and around corners.

\*\*\*This product is available exclusively to Authorized Pride Brick & Mortar Providers. No internet sales are permitted. Available in store only. Other conditions apply.

**400lbs.**  
Weight Capacity

Up To  
**7.3mph**  
Maximum Speed

**43"**  
Turning Radius

Up To  
**21.2**  
Miles Per Charge at 200 lbs.





# NAVIGATE TIGHT CORNERS

## ZERO TURN8

Say goodbye to ordinary, four-wheel scooters! It's time to turn a corner and say hello to stunning innovation and design. The **Zero Turn 8** (*Jazzy® Zero Turn*) is a four-wheel scooter that offers a three-wheel turning radius. Tight corners, narrow aisles and small spaces are no match against the Zero Turn 8, which is loaded with features to make your ride comfortable and stylish.

# ZERO TURN8

Sapphire Blue

Pewter

Pearl White

Black Onyx

Green  
Moonstone

Fire Opal



Our patent-pending, intelligent-turning *iTurn Technology™* ensures effortless navigation in small spaces and around corners.



FDA Class II Medical Device\*

**325 lbs.**

Weight Capacity

Up To  
**5.9 mph**

Maximum Speed

**37.25"**

Turning Radius

Up To  
**13.1**

Miles Per Charge at 200 lbs.







# ELEVATE YOUR LIFE

**Jazzy  
air<sub>2</sub>**

Rise each day and see your world in a whole new way! Empower your life with the sleek and contemporary styling of the Jazzy Air® 2. With up to 12 inches of elevation in just 11 seconds, the world is within reach. Safely drive at 4 mph while elevated. No matter where your day takes you, this power chair gives you the confidence you need!

# Jazzy air<sub>2</sub>

Silver
Pearl White
Ruby Red
Onyx Black
Pink Topaz
Sapphire Blue
Citrine Yellow
Garnet Red (matte)
Tanzanite (matte)
Black Pearl (matte)

Angle and height-adjustable headrest

Angle, height and depth-adjustable armrests

High-back or solid seat pan, reclining sport seat

Storage compartment

Drives safely at 4 mph while elevated at 12"

Integrated foot platform

Bright LED marker lights

Active-Trac<sup>®</sup> Suspension with Mid-Wheel 6<sup>®</sup> Technology

FDA Class II Medical Device\*



**300lbs.**  
Weight Capacity

Up To  
**4mph**  
Maximum Speed  
Elevated and  
Non-Elevated

**24.25"**  
Turning Radius

Up To  
**19.2**  
Miles Per Charge at 200 lbs.







FDA Class II Medical Device\*

# TRAVELER

**Jazzy**<sup>®</sup>  
**carbon**

The Jazzy<sup>®</sup> Carbon uses ultralight but exceptionally strong carbon fiber technology for its frame and a lightweight, reliable lithium-ion battery.

Travelers keep their Jazzy Carbon charged and standing by for spontaneous jet-setting adventures. Weighing in at a stunning 43.6 pounds, the Jazzy Carbon still provides top speeds up to 3.7 mph with a weight capacity of 300 pounds.

**Jazzy**  
carbon

**43.6 lbs.**  
TOTAL WEIGHT  
WITH BATTERIES

Carbon Fiber  
Black

Padded seat back  
with storage pocket

Ultra-lightweight,  
carbon fiber frame

Airline compliant  
lithium battery\*

FDA Class II Medical Device\*



Removable joystick  
with integrated  
USB charger

Compact base easily  
folds and unfolds for  
easy transport in a  
few easy steps

Front suspension  
for a smooth ride

Under-seat  
storage

**300 lbs.**  
Weight Capacity

Up To  
**3.7 mph**  
Maximum Speed

**24"**  
Turning Radius

Up To  
**9.3**  
Miles Per Charge



Height: 28.3" | Width: 24" | Depth: 12.2"







FDA Class II Medical Device\*

# GET OUT THERE!

## GOCHAIR

The night is young and nothing is stopping you. From chic restaurants to fabulous night clubs, the **Go Chair®** is your ticket to unforgettable experiences. This compact power chair is lightweight yet offers the comfort you need for hours of fun. The Go Chair® disassembles easily and offers built-in storage to keep everything you need close at hand. Embrace the nightlife and get out there with the Go Chair®!

- Ruby Red
- Pearl White
- Amber Orange
- Citrine Yellow
- Rose Quartz
- Sapphire Blue



Height and width adjustable armrests

Five-point swivel seat with latch release

Dual swing-away, under seat storage bins

One-hand feather-touch disassembly

Large footplate

FDA Class II Medical Device\*


**300 lbs.**  
Weight Capacity

Up To  
**3.7 mph**  
Maximum Speed

Up To  
**25.5"**  
Turning Radius

Up To  
**13.2**  
Miles Per Charge at 200 lbs.





1. Cup Holder | 2. Lap Belt | 3. Military Patches | 4. Oxygen Holder | 5. Rear Basket | 6. RAM® X-Grip Cell Phone Holder  
 7. Saddlebag | 8. Safety Flag | 9. Walker Holder | 10. Weather Cover | 11. XLR USB Charger

# ACCESSORIES

Scan to check out our full range of accessories.



Scooter/Travel Accessories



Jazzy Power Chair Accessories



MODEL NAME / NUMBER	GO GO® ENDURANCE Li / S54LXLIT	ZERO TURN 10 / S710ZT	ZERO TURN 8 / JAZZYZTS
WEIGHT CAPACITY	325 lbs.	400 lbs.	325 lbs.
MAXIMUM SPEED <sup>1</sup>	Up to 4.2 mph	Up to 7.3 mph	Up to 5.9 mph
GROUND CLEARANCE <sup>3</sup>	1.5" at motor, 1.5" at anti-tip, 3.25" at side deck	2" at motor, 3" at mid-frame	1.3" at motor, 1.5" at anti-tip, 3.5" at side deck
TURNING RADIUS <sup>3</sup>	48.25"	43"	37.25"
OVERALL LENGTH / WIDTH <sup>3</sup>	42" / 20.5"	48" / 24.25"	45.75" / 20.75"
SEAT-TO-GROUND HEIGHT RANGE	21.25"-22.75"	23.25"-24.75"	22"-23.5"
SEAT-TO-DECK HEIGHT RANGE	16.25"-17.75"	17.75"-19.25"	17"-18.5"
FRONT WHEELS (TIRES)	2" x 8" solid	9" solid	2" x 7" solid
REAR WHEELS (TIRES)	2.5" x 8" solid	10.75" solid	2.5" x 9" solid
DRIVE WHEELS	N/A	N/A	N/A
SUSPENSION	Front and rear (CTS Suspension)	Front and rear (CTS Suspension)	Front and rear (CTS Suspension)
RANGE PER CHARGE <sup>1,2</sup>	Up to 7.2 miles (8AH) Up to 14.4 miles at (16 AH MAXX)	Up to 15.9 miles at 400 lbs. Up to 21.2 miles at 200 lbs.	Up to 11 miles at 325 lbs. Up to 13.1 miles at 200 lbs.
HEAVIEST PIECE <sup>3</sup>	40.75 lbs. (front section)	75 lbs. (front section)	48 lbs. (rear section)
STANDARD SEATING	Type: Compact, foldable Weight: 23 lbs. (17"x17") Weight: 27.75 lbs. (20"x18") Material: Black vinyl Dimensions: 17"x17"; 20"x18"	Type: Foldable, deluxe contour Weight: 30 lbs. Material: Black vinyl Dimensions: Width: 18", Depth: 17"	Type: Compact, foldable Weight: 27 lbs. Material: Black vinyl Dimensions: Width: 17", Depth: 17"
DRIVE SYSTEM	24-volt DC motor, sealed transaxle	24-volt DC motor, rear-wheel drive, dual in-line	24-volt DC motor, rear-wheel drive, dual in-line
DRIVE CONTROLLER	N/A	N/A	N/A
BRAKING SYSTEM	Regenerative and electromechanical	Regenerative and electromechanical	Regenerative and electromechanical
BATTERIES <sup>5</sup>	8 AH lithium ion pack or 16 AH MAXX lithium ion pack	(2) 40 AH	(2) 22 AH (AGM)
BATTERY WEIGHT <sup>4</sup>	8 AH – 8 lbs. 16 AH – 12 lbs.	30.25 lbs. each	32 lbs. battery box
BATTERY CHARGER	Off-board, 3.5 amp (lithium-ion)	Off-board, 3.5 amp	Off-board, 3.5 amp
HCPCS CODES	K0899**	K0899**	K0899**
WARRANTY	Frame: lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Battery: 2-year limited	Frame: lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Battery: 6-month-limited	Frame: lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Battery: 6-month-limited

# SPECIFICATIONS

MODEL NAME / NUMBER	JAZZY AIR® 2 / JAIR21001	JAZZY® CARBON	GO CHAIR® / GOCHAIR
WEIGHT CAPACITY	300 lbs.	300 lbs.	300 lbs.
MAXIMUM SPEED <sup>1</sup>	Up to 4 mph	Up to 3.7 mph	Up to 3.7 mph
GROUND CLEARANCE <sup>3</sup>	1"	2"	1.7" at anti-tip, 1.8" at motor, 2.7" at footplate
TURNING RADIUS <sup>3</sup>	24.25"	24"	25.5" with foot platform
OVERALL LENGTH / WIDTH <sup>3</sup>	41.93" / 25.25"	37" / 24"	33" with foot platform / 22.5"
SEAT-TO-GROUND HEIGHT RANGE	21.25" - 33"	20.25" (18.75" without cushion)	22" - 23"
SEAT-TO-DECK HEIGHT RANGE	N/A	N/A	16.25" - 17.25" (platform)
FRONT WHEELS (TIRES)	6" solid	7" solid casters	6" solid casters
REAR WHEELS (TIRES)	6" solid	N/A	9" solid (drive)
DRIVE WHEELS	10" solid	8.5" solid	N/A
SUSPENSION	Active-Trac® with Mid-Wheel 6® Technology	Front	Limited
RANGE PER CHARGE <sup>1,2</sup>	Up to 16.8 miles (standard) Up to 19.2 miles (optional)	Up to 9.3 miles	Up to 8.7 miles (standard) at 300 lbs. Up to 13.2 miles at 200 lbs.
HEAVIEST PIECE <sup>3</sup>	202 lbs. w/seat (no disassembly)	39.6 lbs. (base)	36 lbs. (rear)
STANDARD SEATING	High-back sport seat Material: Black vinyl Dimensions: 16" x 16" - 18"; 18" x 18" - 20"; 20" x 20" - 22"	Material: Black, Foam Dimensions: 17" x 17"	Type: Compact, foldable Material: Black vinyl Dimensions: 18"x17"-30 lbs., 20"x18"-32 lbs. Weight: 18"- 30 lbs., 20"- 32 lbs.
DRIVE SYSTEM	Two-motor, Mid-Wheel 6® drive	24v/250w brushless motors	Two-motor rear-wheel drive
DRIVE CONTROLLER	40A, PG VR2	35A, Brushless Micon	40A, PG GC3
BRAKING SYSTEM	Regenerative and electromechanical	Regenerative and electromechanical	Regenerative and electromechanical
BATTERIES <sup>5</sup>	(2) U1 (standard) (2) 40 AH (optional)	24V 12 AH (288WH) lithium ion	(2) 18 AH
BATTERY WEIGHT <sup>4</sup>	U1 – 24.5 lbs each; 40 AH – 32.5 lbs. each	4 lbs.	14 lbs. each 29.5 lbs. battery box
BATTERY CHARGER	Off-board, 5 amp	2 amp, lithium-ion off-board	2 amp, off-board
HCPCS CODES	K0899**	K0899**	K0899**
WARRANTY	5-year limited on frame Drive motors: 13-month limited Electronics: 13-month limited Battery: 13-month limited	Lifetime limited on frame 1-year limited on electrical control/ lever system, motor/drive system, bearing sleeve and battery	5-year limited on frame Drive motors: 13-month limited Electronics: 13-month limited Battery: 13-month limited

# SPECIFICATIONS

# PRIDE<sup>®</sup>

MOBILITY  
LIVE YOUR BEST<sup>®</sup>

(US) 800-800-8586 | (CAN) 888-570-1113 | [pridemobility.com](http://pridemobility.com) |    



**Note:** Transport of the Jazzy<sup>®</sup> Carbon item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Jazzy<sup>®</sup> Carbon and its lithium-ion battery are permitted.

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
- 2 Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.
- 5 AGM or gel cell type required.

**NOTE:** All specifications subject to change without notice.  
The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.

\*Pride<sup>®</sup> FDA Class II Medical Devices are designed to aid individuals with mobility impairments.

\*\*K0899 – power mobility device, not coded by dme pdac or does not meet criteria.

\*\*\*NM – This product is not a medical device and is not intended to assist, treat, diagnose, or alleviate any medical condition or disability.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

